


# Alaris<sup>®</sup> Syringe Pumps

## Technical Service Manual



***This manual has been prepared for use by qualified service personnel only.  
Cardinal Health cannot accept any liability for any breakdown or deterioration in  
performance of parts or equipment resulting from unauthorised repair or modification.***

**🏢 Cardinal Health, 1180 Rolle, Switzerland**



**Alaris®, Guardrails® and Asena® are registered trademarks of  
Cardinal Health, Inc. or one of its subsidiaries.**

**All other trademarks belong to their respective owners.**

# ***Contents***

## **Chapter**

|                                             |           |
|---------------------------------------------|-----------|
| <b>1. General Information</b>               | <b>4</b>  |
| <b>2. Configuration &amp; Calibration</b>   | <b>9</b>  |
| <b>3. Routine Maintenance</b>               | <b>21</b> |
| <b>4. Troubleshooting</b>                   | <b>30</b> |
| <b>5. Circuit Descriptions</b>              | <b>37</b> |
| <b>6. Spare Part Replacement Procedures</b> | <b>40</b> |

## **Appendix**

|                                                     |           |
|-----------------------------------------------------|-----------|
| <b>A. Specifications</b>                            | <b>62</b> |
| <b>B. Disposal</b>                                  | <b>69</b> |
| <b>C. Spare Parts Listing</b>                       | <b>71</b> |
| <b>D. Fitting &amp; Replacement Guidelines</b>      | <b>75</b> |
| <b>E. Configuration &amp; Drug Protocol Records</b> | <b>78</b> |
| <b>F. Service Contacts</b>                          | <b>85</b> |
| <b>G. Document History</b>                          | <b>87</b> |

# ***Chapter 1***

## ***General Information***

### **In this chapter**

|                                     |          |
|-------------------------------------|----------|
| <b>Introduction</b>                 | <b>5</b> |
| <b>General precautions</b>          | <b>5</b> |
| <b>Front panel and main display</b> | <b>6</b> |
| <b>Controls and indicators</b>      | <b>6</b> |
| <b>Loading a syringe</b>            | <b>7</b> |
| <b>Starting the pump</b>            | <b>8</b> |
| <b>Basic features</b>               | <b>8</b> |

## Introduction

The Alaris® Syringe Pumps are designed to deliver a continuous and accurate infusion whenever small fluid volumes need to be administered with great precision. High performance, comprehensive alarm protection and sophisticated monitoring systems combined with simple operation make these syringe pumps suitable for both general and critical infusions in a variety of areas within a hospital.

The Asena® brand name has been recently changed to the Alaris® brand name. This change in brand name has no effect on the intended use or functionality of the product. Recommended disposable products for use with this product may refer to either the Asena® brand name or Alaris® brand name and both types are suitable for use with this infusion pump.

### Familiarity

Ensure that you are fully familiar with this syringe pump by carefully studying the Directions for Use (DFU) prior to attempting any repairs or servicing.

As part of a policy of continuous improvement, product enhancements and changes are introduced from time to time.

### Purpose

This Technical Service Manual shows how to set up, test and maintain the following Alaris® Syringe Pump models:

- Alaris® GH Syringe Pump
- Alaris® GH Syringe Pump with Guardrails® Safety Software
- Alaris® CC Syringe Pump
- Alaris® CC Syringe Pump with Guardrails® Safety Software
- Alaris® PK Syringe Pump
- Alaris® TIVA Syringe Pump
- Alaris® GS Syringe Pump

It is intended for use by personnel experienced in medical equipment testing and maintenance procedures.

### Symbols



Wherever you see this symbol in the manual you will find a Hints & Tips note that we hope you will find useful. These notes provide useful advice or information that may help you perform the task more effectively.



Wherever you see this symbol in the manual you will find a Toolbox note that highlights an aspect of test or maintenance that is important for you to know about. A typical example is a software upgrade that you should check has been installed.

## General precautions



Please read the general Operating Precautions described in the Directions for Use carefully prior to using the pump.



This pump contains static-sensitive components. Observe strict precautions for the protection of static sensitive components when attempting to repair and service the pump.




An explosion hazard exists if the pump is used in the presence of flammable anaesthetics. Exercise care to locate the pump away from any such hazardous sources.



An electrical shock hazard exists if the pumps casing is opened or removed. Refer all servicing to qualified service personnel.



This pump is protected against the effects of high energy radio frequency emissions and is designed to fail safe if extremely high levels of interference are encountered. Should false alarm conditions be encountered, either remove the source of the interference or regulate the infusion by another appropriate means.




If the pump is dropped, subjected to excessive moisture, humidity or high temperature, or otherwise suspected to have been damaged, remove it from service for inspection by qualified service personnel.



When connected to an external power source, a three-wire (Live, Neutral, Earth) supply must be used. If the integrity of the external protective conductor in the installation or its arrangement is in doubt, the pump should be operated from the battery.

## Front panel and main display

The display shown is for general guidance only. For pump specific front panel and main display information refer to relevant Directions For Use.



\* "Profile" is only available on an Alaris® Syringe Pump with Guardrails® Safety Software with a Data Set loaded.

\*\* Pressure Information is only displayed on the Alaris® CC Syringe Pump.

## Controls and indicators



**ON/OFF button** - Press once to switch the pump ON. Press and hold down for 3 seconds to switch the pump OFF.



**RUN button** - Press to start the infusion. The Green LED will flash during infusion.



**HOLD button** - Press to put the infusion on hold. The amber LED will be lit while on hold.



**MUTE button** - Press to silence alarms.



**PURGE/BOLUS button** - Press to access **PURGE** or **BOLUS** soft keys. Press and hold down soft key to operate.

**PURGE** the extension set during set up.

- Pump is on hold
- Extension set is not connected to the patient
- Volume Infused (VI) is not added

**BOLUS** fluid or drug delivered at an accelerated rate.

- Pump is infusing
- Extension set is connected to the patient
- VI is added



**OPTION button** - Press to access optional features.



**PRESSURE button** - Press to display the pumping pressure and alarm level.




**BLANK SOFTKEYS** - Use in conjunction with the prompts shown on the display.



**CHEVRON keys** - Double or single for faster/slower, increase or decrease of values shown on main display.



**BATTERY indicator** - When illuminated, indicates that the pump is running on the internal battery. When flashing, indicates that the battery power is low, with less than 30 minutes of use remaining.




**AC POWER indicator** - When illuminated, indicates that the pump is connected to an AC power supply and the battery is being charged.

## Loading a syringe

1. Squeeze the finger grips together on the plunger holder and slide the mechanism to the right. Pull the syringe clamp forward and down.
2. Place the syringe barrel flange in the slot between the two blue sections of the syringe flange clamp.





**Note:** Earlier pumps have a single section syringe flange clamp. In this instance place the syringe barrel flange in the space between the syringe clamp and the syringe flange clamp.





**To ensure the syringe is loaded correctly, check the syringe remains in position before the syringe clamp is closed.**

3. Lift the syringe clamp until it locks against the syringe barrel.
4. Squeeze the finger grips on the plunger holder and slide the mechanism to the left until it reaches the plunger end.
5. Release the finger grips. Ensure that the plunger grippers are securing the plunger in place and the upper finger grip returns to its original position.



**Ensure that both plunger grippers are fully locked onto the plunger flange and the upper finger grip has returned to its original position.**

## Guide to handling the pressure disc



**When withdrawing the pressure disc, make sure you pull it (with your finger inside the disc recess) towards the front of the pump as shown in the diagram above.**

## Starting the pump

1. Push pump on to bar, the Alaris® Docking Station or mount to pole. Connect to AC Mains.
2. Press the  button to switch pump ON.
3. Follow SETUP / Drug / Profile instructions as per Directions For Use.
4. Load syringe.
5. Insert pressure disc into pressure transducer (Model CC only, required for dedicated mode).
6. Confirm syringe.
7. Change the rate if necessary using the   keys.
8. Purge: Press the  button followed by the **PURGE** softkey.
9. Connect the pump to test equipment as required (see chapter 2 Configuration & Calibration).
10. Press the  button to start the infusion.

## Basic features

|                                                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|--------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Pressure Level</b><br>(not available on Model GS)               | Press the  button. Pressure Alarm level and current pressure level are shown on graph. Use the   keys to adjust pressure alarm level.                                                                              |
| <b>Volume to be Infused</b><br>(not available on Models TIVA & GS) | Press <b>VTBI</b> . Enter VTBI using the   keys. Press <b>OK</b> . Select rate at end of VTBI using the   keys. Press <b>OK</b> . |
| <b>Clear Volume Infused</b>                                        | Press <b>VOLUME</b> . Select <b>YES</b> (to clear) or <b>NO</b> (to exit).                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Purge/Bolus</b>                                                 | Press the  button, use the   keys to set rate, then press and hold down the <b>PURGE</b> or <b>BOLUS</b> softkey.                                                                                                  |
| <b>Hands-free Bolus*</b><br>(Model TIVA only)                      | Press the  button, if required, select <b>HANDSFREE BOLUS</b> , use the   keys to adjust bolus dose. Use rate softkey to adjust bolus rate if required. Press <b>BOLUS</b> softkey to deliver.                 |



**During PURGE/BOLUS, the pressure limit alarms are temporarily increased to maximum levels.**

## Option Button

|                           |                                                                                                          |
|---------------------------|----------------------------------------------------------------------------------------------------------|
| <b>Drug Library*</b>      | Select Drugs and Dosing (Models CC & TIVA) or Drug name (Model GH).                                      |
| <b>Set VTBI Over Time</b> | Specify a VTBI and delivery time ( <i>not available on Models TIVA &amp; GS</i> ). Pump must be on hold. |
| <b>24 Hour Log</b>        | Volume infused log over 24 hours and accumulative total (on the hour only).                              |
| <b>Rate Lock</b>          | Disables changing rate once infusion has started ( <i>not available on Model TIVA</i> ).                 |



**When rate lock is enabled, the following are unavailable:**

1. Rate changes/titration
2. Bolus/Purge
3. Switch pump off

|                       |                                                                                                                                        |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Event Log*</b>     | Access Event Log. The Event Log holds up to 1500 individual events. Pumps with Guardrails® Software enabled retain one year of events. |
| <b>Hospital Name*</b> | Displays the name of Hospital/Ward/Department as set up on the pump. Accessible while the pump is infusing.                            |

Additional options may be shown, please refer to the relevant Directions For Use, for more information.

\* **Note:** For Guardrails® Software enabled pumps or the Alaris® PK Syringe Pump these options may vary or will not be available. Please refer to the relevant Alaris® Syringe Pump Directions For Use or Guardrails® Editor Directions For Use for comprehensive information.

## ***Configuration & Calibration***

### **In this chapter**

|                                     |           |
|-------------------------------------|-----------|
| <b>Access codes</b>                 | <b>10</b> |
| <b>Dedication options (301/302)</b> | <b>10</b> |
| <b>Data Set Activation (612)</b>    | <b>10</b> |
| <b>Handsfree Bolus (175)</b>        | <b>10</b> |
| <b>Power Lock (711)</b>             | <b>10</b> |
| <b>Configuration options (251)</b>  | <b>11</b> |
| <b>Calibration procedures (243)</b> | <b>15</b> |

## Access codes

The syringe pump software contains a number of configuration and test routines that can be accessed by the user. The majority of tests are 'MENU' driven from a technical access code (see below).

| Code | Description                                                                                                  |
|------|--------------------------------------------------------------------------------------------------------------|
| 123  | Self Test Procedure                                                                                          |
| 166  | External Reprogramming                                                                                       |
| 167  | Teach Learn Procedure                                                                                        |
| 175  | Handsfree Bolus                                                                                              |
| 243  | Calibration Selection Menu                                                                                   |
| 251  | User Configuration Menu                                                                                      |
| 301  | Fully Dedicated                                                                                              |
| 302  | Semi-dedicated                                                                                               |
| 376  | Service Access Menu                                                                                          |
| 401  | Upload Data Set to Pump (Guardrails® Software enabled Pumps and the Alaris® PK Syringe Pump)                 |
| 402  | Download CQI Event Log from Pump (Guardrails® Software enabled Pumps only)                                   |
| 418  | Alternative Alarm Tone. (Not available for Guardrails® Software enabled Pumps & the Alaris® PK Syringe Pump) |
| 499  | Download Data Set from Pump (Guardrails® Software enabled Pumps and the Alaris® PK Syringe Pump)             |
| 611  | Cold Start (RAM Clear)                                                                                       |
| 612  | Data Set activation (Alaris® PK Syringe Pump)                                                                |
| 711  | Power Lock (Alaris® PK Syringe Pump)                                                                         |

Each MENU (and some unique items) has its own three-digit access code that can be entered using the following procedure:

1. Hold down  and turn the pump ON.
2. Enter the required access code using the   keys and the **NEXT** softkey.
3. When the required code shows on screen, press **OK** to confirm.

## Dedication options (301/302)

Fully Dedicated (set using access code **301**) will remind a user that a pressure disc must be fitted to start any infusion. In this mode occlusion pressures are always displayed in mmHg.

Semi-Dedicated (set using access code **302**) will remind a user that a pressure disc must be fitted when drugs and dosing features are used. When a pressure disc is not in use, pressure levels L-0 to L-10 will be displayed.

## Data Set Activation (612)

This code is used to load the predefined pump configuration and drug setup into the non-volatile storage. It is necessary to enter the code **612** after a cold start (code **611**); the configuration and drug setup will then be available in normal operation.

Alternatively a data set may be uploaded as appropriate. See directions for use contained within the Alaris® PK Editor Software package.

## Handsfree Bolus (175)

Enable or disable the Handsfree Bolus. If enabled pressing bolus button displays screen prompting for hands free or hands on. Default volume after clear setup is 0.0. Upper amount restricted to bolus volume limit in general options or drug protocol bolus volume limit.

## Power Lock (711)

Available on the Alaris® PK Syringe Pumps with software V2.3.11 & above.

|                 |                                                                                                                                                                                                                                     |
|-----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Disabled</b> | The new alternative Power Down sequence now allows the user to Power Down the pump whilst the infusion is suspended (on hold) in TCI mode and predictive TIVA mode.                                                                 |
| <b>Enabled</b>  | The Power Down sequence (Power Lock) remains the same where the user may only Power Down the pump by stopping the infusion, selecting 'new operation' from the options menu, confirming the selection, then Powering Down the pump. |

## Configuration options (251)

Enter access code **251** to display the User Configuration menu:

- Drug Library\*** Set drug names list on a Model GH - Select Character Group   Select Character    
To go to next Character use **NEXT**.  
Set drug names and protocols for Models CC and TIVA (see *drug protocol setup instructions on following pages*).
- General Options\*** See table on the following page.
- Clock Set** Set the current date and time. To set the clock, use   and **NEXT** to adjust and **OK** to store.
- Hospital Name\*** Enables establishment name (max 20 characters) to be displayed during the power-up sequence.  
To set the hospital name, use   and **NEXT** to adjust and **OK** to store.
- Enable Syringes\*** Configure the type and size of syringes permitted for use.  
To enable syringes, use   and **SELECT**, to enable/disable and **OK** to store.
- Language** Configure the language used for messages shown on display.  
Select language required using   and **SELECT** to store.
- Contrast** Set the display panel contrast. Use   to adjust contrast and **OK** to store.
- Enable Units\*** Select the type of units permitted for use on the pump. To enable units, use   and **MODIFY**, to enable/disable and **OK** to store.

\* **Note:** For Guardrails® Software enabled pumps and the Alaris® PK Syringe Pump these options may vary or will not be available. Please refer to the relevant Alaris® Syringe Pump Directions For Use or Guardrails® Editor Directions For Use for comprehensive information.

## Configuration options (251) (continued)

### Alaris® TIVA Syringe Pump drug protocol setup

Select Drug Library from Configuration Options (251).

Use   to select drug and press **MODIFY** to modify selected drug or **NEW** to create new drug name.

**QUIT** will return to **251** main menu.

When modifying a drug protocol, pressing **BACK** at any time will take you to the previous step.

#### Modify - Existing drug

**ENABLE/DISABLE** - Enables or disables the drug being available.

**DELETE** - Select Yes to delete from drug library.

**EDIT** - See table below.

#### Edit Drug Protocol - New or existing drug

Press **OK** softkey to confirm each step.

| Drug Option           | To Adjust<br>(Softkeys are shown in <b>Bold</b> )                                                                                                                                                                                                                                                                                                                                                                                         |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Drug Name             | Select Character Group  <br>Select Character  <br>To go to next Character <b>NEXT</b> |
| Concentration Units   |                                                                                                                                                                                                                                                                     |
| Minimum Concentration |   or <b>OFF</b>                                                                                                                                                                                                                                                     |
| Default Concentration |   or <b>OFF</b>                                                                                                                                                                                                                                                     |
| Maximum Concentration |   or <b>OFF</b>                                                                                                                                                                                                                                                     |
| Dose Rate Units       |                                                                                                                                                                                                                                                                     |
| Induction Dose        |   or <b>OFF</b>                                                                                                                                                                                                                                                     |
| Induction Time        |                                                                                                                                                                                                                                                                     |
| Pause After Induction | <b>MODIFY</b>                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Maintenance Rate      |                                                                                                                                                                                                                                                                     |
| Bolus Dose            |                                                                                                                                                                                                                                                                     |
| Bolus Rate            | <b>RATE</b>                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Hands Free Bolus      | <b>MODIFY</b>                                                                                                                                                                                                                                                                                                                                                                                                                             |

## Configuration options (251) (continued)

### Alaris® CC Syringe Pump\* drug protocol setup

Select Drug Library from Configuration Options (251).

Use   to select drug and press **MODIFY** to modify selected drug or **NEW** to create new drug name.

**QUIT** will return to **251** Main menu.

When modifying a drug protocol, pressing **BACK** at any time will take you to the previous step.

#### Modify - Existing drug

**ENABLE/DISABLE** - Enables or disables the drug being available.

**DELETE** - Select Yes to delete from drug library.

**EDIT** - See table below.

#### Edit Drug Protocol - New or existing drug

Press **OK** softkey to confirm each step.

\* **Note:** For Guardrails® Software enabled pumps this option will not be available. Please refer to the relevant Alaris® Syringe Pump Directions For Use or Guardrails® Editor Directions For Use for comprehensive information.

| Drug Option           | To Adjust<br>(Softkeys are shown in <b>Bold</b> )                                                                                                                                                                                                                                                                                                                                                                                         |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Drug Name             | Select Character Group  <br>Select Character  <br>To go to next Character <b>NEXT</b> |
| Dose Rate Units       |                                                                                                                                                                                                                                                                     |
| Maximum Dose          |   or <b>OFF</b>                                                                                                                                                                                                                                                     |
| Default Dose          |   or <b>OFF</b>                                                                                                                                                                                                                                                     |
| Minimum Dose          |   or <b>OFF</b>                                                                                                                                                                                                                                                     |
| Concentration Units   |                                                                                                                                                                                                                                                                     |
| Minimum Concentration |   or <b>OFF</b>                                                                                                                                                                                                                                                     |
| Default Concentration |   or <b>OFF</b>                                                                                                                                                                                                                                                     |
| Maximum Concentration |   or <b>OFF</b>                                                                                                                                                                                                                                                     |
| Maximum Bolus         |   or <b>OFF</b>                                                                                                                                                                                                                                                     |
| Bolus Rate            |                                                                                                                                                                                                                                                                     |
| Pressure Alarm        |   or <b>OFF</b>                                                                                                                                                                                                                                                     |

## Configuration options (251) (continued)

### General options

| Option               | Description                                                                                                  | Models |     |     |      |
|----------------------|--------------------------------------------------------------------------------------------------------------|--------|-----|-----|------|
|                      |                                                                                                              | GS     | GH* | CC* | TIVA |
| NURSE CALL FITTED    | Enables Nurse Call (hardware option).                                                                        | ✓      | ✓   | ✓   | ✓    |
| NURSE CALL INVERT    | When enabled, the nurse call output is inverted.                                                             | ✓      | ✓   | ✓   | ✓    |
| RS232 SELECTED       | Sets the pump's communications to use RS232 (hardware option).                                               | ✓      | ✓   | ✓   | ✓    |
| NEOI WARNING         | Sets the Near End Of Infusion (NEOI) warning time.                                                           | ✓      | ✓   | ✓   | ✓    |
| EOI POINT            | Sets the End Of Infusion volume.                                                                             | ✓      | ✓   | ✓   | ✓    |
| KVO AT EOI           | Enables pump to run at the Keep Vein Open (KVO) rate when End Of Infusion (EOI) is reached.                  | ✓      | ✓   | ✓   | ✓    |
| KVO RATE             | Sets the Keep Vein Open (KVO) rate.                                                                          | ✓      | ✓   | ✓   | ✓    |
| BACK OFF             | Motor will reverse to relieve line pressure when an occlusion occurs.                                        | ✓      | ✓   | ✓   | ✓    |
| AUTO SAVE            | When disabled, the patient information is cleared on power up.                                               | ✓      | ✓   | ✓   | ✗    |
| RATE LOCK            | When enabled, the rate can be locked to prevent unwanted changes of the set infusion rate.                   | ✓      | ✓   | ✓   | ✗    |
| QUIET MODE           | When enabled, the button beeps are muted.                                                                    | ✓      | ✓   | ✓   | ✗    |
| AC FAIL              | When enabled, the AC Power Failure Alarm will sound if AC power is disconnected.                             | ✓      | ✓   | ✓   | ✓    |
| RATE TITRATION       | When enabled, the rate can be changed whilst the pump is infusing.                                           | ✗      | ✓   | ✓   | ✗    |
| PRESSURE DISPLAY     | Enables / disables the Pressure Icon on the main display.                                                    | ✓      | ✓   | ✓   | ✓    |
| AUTO PRESSURE        | Enables / disables the automatic pressure alarm level option.                                                | ✗      | ✗   | ✓   | ✗    |
| AUTO SET PRESSURE    | Automatically sets the line occlusion pressure alarm level to a specified amount above the current pressure. | ✗      | ✗   | ✓   | ✗    |
| AUTO OFFSET          | Adjusts the automatic offset value used by auto pressure and auto set pressure.                              | ✗      | ✗   | ✓   | ✗    |
| HANDS FREE BOLUS     | Enables / disables hands-free bolus.                                                                         | ✗      | ✗   | ✗   | ✓    |
| CAP PRESSURE         | Sets the maximum pressure limit.                                                                             | ✗      | ✓   | ✗   | ✗    |
| PRESSURE DEFAULT     | Sets the default occlusion alarm level.                                                                      | ✓      | ✓   | ✓   | ✓    |
| DEFAULT BOLUS VOLUME | Sets the default hands-free bolus volume for No Drug mode only.                                              | ✗      | ✗   | ✗   | ✓    |
| MAX PRESSURE         | Sets the maximum pressure limit.                                                                             | ✗      | ✗   | ✓   | ✗    |
| WEIGHT               | Sets the default patient weight in kg.                                                                       | ✗      | ✗   | ✓   | ✓    |
| CAP RATE             | Sets the maximum value for infusion rate.                                                                    | ✓      | ✓   | ✓   | ✗    |
| PURGE RATE           | Sets the purge rate.                                                                                         | ✓      | ✓   | ✓   | ✓    |
| PURGE VOLUME LIMIT   | Sets the maximum permissible purge volume.                                                                   | ✓      | ✓   | ✓   | ✓    |
| PURGE SYRINGE        | Prompt to purge syringe after confirmation.                                                                  | ✓      | ✓   | ✓   | ✓    |
| BOLUS                | Enables / disables the bolus feature.                                                                        | ✓      | ✓   | ✓   | ✗    |
| DEFAULT BOLUS        | Sets the default bolus rate.                                                                                 | ✓      | ✓   | ✓   | ✓    |
| CAP BOLUS RATE       | Sets the maximum value for bolus rate.                                                                       | ✓      | ✓   | ✓   | ✗    |
| BOLUS VOL LIMIT      | Sets the maximum permissible bolus volume.                                                                   | ✓      | ✓   | ✓   | ✗    |
| MANUAL BOLUS         | Volume infused will be increased if plunger is manually moved in and syringe remains confirmed.              | ✓      | ✓   | ✓   | ✓    |
| CALL BACK TIME       | Adjusts the time for the pump to sound the call back alarm.                                                  | ✓      | ✓   | ✓   | ✓    |
| VTBI CLEAR RATE      | Rate will be set to zero when VTBI has been set-up with stop as the end rate.                                | ✗      | ✓   | ✓   | ✗    |
| EVENT LOG DISPLAY    | Enables / disables the event log display.                                                                    | ✓      | ✓   | ✓   | ✓    |
| BATTERY ICON         | Enables / Disable the Battery Icon on the main display.**                                                    | ✓      | ✓   | ✓   | ✓    |
| AUDIO VOLUME         | Sets the alarm volume of the pump at high, medium or low.                                                    | ✓      | ✓   | ✓   | ✓    |
| AUTO NIGHT MODE      | Sets Backlight to dim between 21:00 and 06:00hrs.                                                            | ✓      | ✓   | ✓   | ✓    |

#### Key:

- ✓ = available option
- ✗ = unavailable option

\* For Guardrails® Software enabled pumps and the Alaris® PK Syringe Pump these options may vary or will not be available, with only the first three options listed in table above adjustable in the General Options on the pump. Please refer to the relevant Alaris® Syringe Pump Directions For Use or Guardrails® Editor Directions For Use for comprehensive information.


\*\* For Alaris® GS Syringe Pump the battery icon can be seen via the Options menu by pressing the  key.

## Calibration procedures (243)

Enter access code **243** to display the Calibration Selection menu (see Access Codes).

### SYRINGE CLAMP calibration

Fit calibration tool into position on pump as shown in Steps 1-2 and close the clamp. At each step, CAL is displayed if value is within tolerances. Press CAL button to store calibration point.




Calibration tool required: 1000TG00095


**Note:** If CAL is not displayed, check for correct positioning of calibration tool. If calibration cannot be performed, repairs to pump may be necessary.

**Note:** The calibration values shown on the displays are for illustrative use only and may vary.


#### Step 1



#### Step 2



#### Step 3




## Calibration procedures (243) (continued)

### PLUNGER POS (position) calibration

Fit calibration tool in position on pump as shown in Steps 1-3. At each step CAL is displayed if value is within tolerance. Press CAL button to store calibration point.


**Note:** If CAL is not displayed, check for correct positioning of calibration tool. If calibration cannot be performed, repairs to pump may be necessary.

**Note:** The calibration values shown on the displays are for illustrative use only and may vary.





Calibration tool required: 1000TG00095

#### Step 1




#### Step 2


Close-up of calibration tool, showing locking clip in position.



#### Step 3



#### Step 4




## Calibration procedures (243) (continued)

### SYRINGE FORCE calibration

#### Precondition:

This precondition the mechanism and should only be done if motorplate or chassis has been replaced. Fit Calibration tool as shown, zero the gauge, run until gauge registers 10kgf and then carefully declutch mechanism and withdraw plunger. Do not press any button during this procedure.



Calibration tool required:  
0000TG00200 (top) or  
0000TG00020 (bottom)



**To convert Kilograms of Force (kgf) to Newtons (N) multiply by 9.806650. For example 10kgf = 98.07N.**




**Excessive force will damage the plunger mechanism. Do not apply more than 10 kgf ±0.05kgf to the plunger mechanism at any time.**

**Note:** The calibration values shown on the displays are for illustrative use only and may vary.

#### Precondition

**10kgf ±0.05kgf**




Fit Calibration tool and position plunger as shown in Steps 1 to 3, zero the gauge. At each step press CAL when required calibration force is reached.

**Note:** If CAL does not appear in display, check for correct positioning of tool. If calibration cannot be performed, repairs to pump may be necessary.

Allow 30 seconds for pressure to stabilise following any preconditioning calibration.

#### Step 1

**0kgf ±0.05kgf**




## Calibration procedures (243) (continued)

### SYRINGE FORCE calibration (continued)


#### Step 2

3Kgf  $\pm 0.05$ Kgf




#### Step 3

10Kgf  $\pm 0.05$ Kgf



#### Step 4



**Use of the 0000TG00200 Digital Occlusion Testgear.**

**The 0000TG00200 Occlusion testgear uses a digital force gauge to register applied forces. Please refer to the MecMesin Compact Gauge Operation Instructions supplied for detailed operational information and power options and requirements.**

**To prepare the testgear for use, load into the syringe pump.**

- **Ensure there is nothing touching the testgear plunger (such as the syringe plunger drive).**
- **Turn on the Compact Gauge using the 'On/Zero' key.**
- **Select 'kg' force units, and 'MAX' reading option.**
- **If the display indicates other than 0.00kg, zero the system using the 'On/Zero' key.**

**Operate the system as required for performing the calibration activity.**

**Before the next use, ensure the 'MAX' reading is cleared using the 'On/Zero' key.**

## Calibration procedures (243) (continued)


### LINE PRESSURE calibration – Alaris® CC Syringe Pump ONLY

#### Tools required:

Pressure gauge (range 0-1400 mmHg)  
(Tolerance +/- 0.1mmHg)

Dedicated pressure sensing disc extension set  
(i.e. G30402)


50ml Luer-lock syringe




Load pressure disc infusion set into transducer. Connect infusion set to syringe and gauge. Using syringe, apply pressure required as shown at steps 1-3. At each step press CAL when required calibration pressure is displayed on pressure gauge.

**Note:** The calibration values shown on the displays are for illustrative use only and may vary.


#### Step 1 25mmHg ± 1mmHg




#### Step 2 500mmHg ± 1mmHg



#### Step 3 1000mmHg ± 1mmHg



#### Step 4



## Calibration procedures (243) (continued)

### BATTERY calibration

1. Connect the Pump to AC mains.
2. Select BATTERY CALIBRATION from menu and press **OK**.
3. The pump will automatically run the battery calibration. Battery calibration cycles the battery through a charge, discharge and re-charge sequence during which the gas gauge within the battery pack will be updated with a measurement of the current capacity of the cells.




**Battery compartment should be ventilated during calibration (open battery cover). Pump may fail calibration if too hot, so care should be taken not to calibrate too many pumps in close proximity (in a docking station, for example). Ensure that the battery is supported as you open the battery compartment.**



**Disconnecting the AC mains at any time during calibration will cause battery calibration to fail.**

4. When calibration is complete, the following is shown on the display:



|   | Value                        | Description                                                                                                             | Pass Criteria                         |
|---|------------------------------|-------------------------------------------------------------------------------------------------------------------------|---------------------------------------|
| 1 | Battery Capacity             | Pack capacity value updated after measured discharge phase (if changed).                                                | Greater than 2100mAh                  |
| 2 | Current Battery Charge Level | Current charge in pack.                                                                                                 | n/a                                   |
| 3 | Battery Voltage              | Current pack voltage.                                                                                                   | n/a                                   |
| 4 | Initial Charge Time          | Time taken during initial charge phase. Initial charge phase checks pack is fully charged and if not it is charged.     | Lower than 2 hours 59 minutes         |
| 5 | Discharge Time               | Time taken during measured discharge phase. Pack is discharged to determine how much charge is available from the pack. | Between 4 hours 15 minutes & 10 hours |
| 6 | Final Charge Time            | Time taken during final charge phase. Pack is fully recharged ready for use.                                            | Lower than 2 hours 59 minutes         |

5. All pass criteria (see table above) should be met and the pump should display FINISH at the end of the calibration otherwise calibration has failed. If calibration has failed retry calibration or replace battery.
6. Press **OK** to exit.



**The plunger drive will move automatically during the discharge phase, so ensure that the plunger drive is not obstructed during calibration (remove syringes etc).**

## ***Routine Maintenance***

### **In this chapter**

|                                            |           |
|--------------------------------------------|-----------|
| <b>Self-test procedure (123)</b>           | <b>22</b> |
| <b>Comms test (123)</b>                    | <b>23</b> |
| <b>Data transfer</b>                       | <b>24</b> |
| <b>Information logs (376)</b>              | <b>27</b> |
| <b>Recommended cleaning and inspection</b> | <b>28</b> |
| <b>Performance verification procedure</b>  | <b>29</b> |

## Routine Maintenance

For routine maintenance, the following self-test and performance verification procedures should be performed. Refer to the relevant Directions for Use for the recommended routine maintenance period.

### Self-test procedure (123)

#### Self-tests included in full test

Enter access code **123** to view the Test Selection menu (see Access Codes in chapter 2). Refer to table below for the tests in each menu item.

| Test Section     | Test             | Action                                                                                                                                                                                          |
|------------------|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Software         | Software info    | Displays the software version.                                                                                                                                                                  |
|                  | Data Set Info    | Displays the Data Set information. (pumps with Guardrails® Safety Software only)                                                                                                                |
| Safety Processor | Safety ID        | Check displays the version of the safety ID.                                                                                                                                                    |
|                  | Safety LED       | Check red LED illuminated.                                                                                                                                                                      |
|                  | Safety Alarm     | Check Backup alarm sounds.                                                                                                                                                                      |
| Full only        | Serial Number    | Check displays serial number of unit.                                                                                                                                                           |
|                  | Language         | Check displays correct language.                                                                                                                                                                |
|                  | Real-time Clock  | Check displays correct date and time.                                                                                                                                                           |
|                  | Service Date     | Check displays date when service is next required.                                                                                                                                              |
| Sensor           | Disc Detect      | Check the display changes correctly to indicate if a disc is Out or In (Model CC only).                                                                                                         |
|                  | Line Pressure    | Check pressure is 000mmHg +/-20mmHg with no pressure applied (Model CC only).                                                                                                                   |
|                  | Motor Encoder    | Check motor runs and Passed is displayed.                                                                                                                                                       |
|                  | Drive Engage     | Check display indicates Drive Engaged or Disengaged when clutched/ declutched.                                                                                                                  |
|                  | Plunger Fit      | Check display indicates if the Plunger button is Out or In.                                                                                                                                     |
|                  | Plunger Position | Check display smoothly and continuously changes during full plunger travel.                                                                                                                     |
|                  | Syringe Clamp    | Insert the syringe size calibration tool (1000TG00095) and check the following values are displayed for diameters inserted:<br>12mm diameter = 11.5 to 12.5mm<br>32mm diameter = 31.5 to 32.5mm |
|                  | Syringe Force    | Check motor runs and syringe force is displayed.                                                                                                                                                |
| Battery          | Battery          | Check displays values in CAP, CHR and VOL; no dashes should be seen.                                                                                                                            |
| Audio            | Audio Speaker    | Check the main audible alarm sounds.                                                                                                                                                            |
| Visual Indicator | Display          | Check that all of the display pixels are illuminated.                                                                                                                                           |
|                  | Backlight        | Check that the backlight switches from LOW to HIGH when indicated.                                                                                                                              |
|                  | Battery LED      | Check the Battery LED (Amber) flashes.                                                                                                                                                          |
|                  | Start LED        | Check the Start LED (Green) flashes.                                                                                                                                                            |
|                  | Stop LED         | Check the Stop LED (Amber) flashes.                                                                                                                                                             |
|                  | Warning LED      | Check the Warning LED (Amber) flashes.                                                                                                                                                          |
|                  | Alarm LED        | Check the Alarm LED (Red) flashes.                                                                                                                                                              |
| Key              | Keypad           | Press the key indicated and check changes to next key.                                                                                                                                          |
| Comms            | Comms            | RS232 only. Check Nurse call and RS232 operation.                                                                                                                                               |

## Self-test procedure (123) (continued)

### Self-tests not included in full test

| Test Section        | Test             | Action                                                                        |
|---------------------|------------------|-------------------------------------------------------------------------------|
| Remote              | Remote           | Check the function of the IrDA output for remote access                       |
| Calibration records | Syringe clamp    | Displays calibration values for Closed and Open positions.                    |
|                     | Plunger position | Displays calibration values for Left, Middle and Right positions.             |
|                     | Syringe force    | Displays calibration values for 0, 3 and 10 kgf.                              |
|                     | Line pressure    | Displays calibration values for 25, 500 and 1000mmHg (Model CC only)          |
| Linearity           | Linearity        | Check the mechanism runs full travel and graph displays smooth linear travel. |
|                     | Occlusion base   | Check the occlusion base level is within tolerance shown on graph.            |

## Comms test (123)

Select COMMS TEST from the displayed menu. **Note:** Section only applicable if RS232 Hardware option is fitted.

No specific customer test equipment is available to conduct the RS232 on nurse call alarm tests. It is assumed that the customer will have associated systems that make use of the RS232 and nurse call options, hence:

The nurse call system can be tested, once connected to the customer facility, by running the pump and simulating an alarm condition (e.g. Disengaging the drive while running).

The RS232 system can be tested by communicating with the pump using a customer application.

If no customer systems are available for in-use testing, the following connections to the 9 pin D type output socket will permit testing. It is recommended that all test connections are made via a 9 way D type plug that can be fitted into the pump socket.

#### Test

RS232 Test


#### Description

Only available when Nurse Call option is fitted. **Note:** NURSE CALL FITTED & RS232 SELECTED must be enabled (✓) in access code **251** General Options for this test to work. Connect the 9-pin D type connector to the 9 pin D type output socket at the rear of the pump. The display '\_\_\_\_' will change to PASS if the communications test is successful.

Nurse Call

Only available when Nurse Call option is fitted. **Note:** NURSE CALL FITTED & RS232 SELECTED must be enabled (✓) in access code **251** General Options for this test to work. Locate the 9-pin D type connector at the rear of the pump. Check that the Nurse Call relay switches from NC to NO connections during test.

RS232 pinout



| Pin Number<br>(Pump Socket Id) | Required Action                               | Comments                                                                                             |
|--------------------------------|-----------------------------------------------|------------------------------------------------------------------------------------------------------|
| 1                              | Nurse call relay - normally closed connection | With nurse call test in progress - Confirm continuity with pin 5 - Alternately switches with pin 8.  |
| 2                              | Link pin 2 to pin 3                           | RS232 Tx and RX test link. With RS232 test in progress - Confirm 'PASS' is displayed on test screen. |
| 3                              | See pin 2                                     | -----                                                                                                |
| 4                              | Not used                                      | -----                                                                                                |
| 5                              | 0 volt line                                   | With respect to pin 7.                                                                               |
| 6                              | Not used                                      | -----                                                                                                |
| 7                              | Apply 9 to 12 volts DC                        | RS232 Power supply - with respect to pin 9.                                                          |
| 8                              | Nurse call relay - normally open connection   | With nurse call test in progress - Check continuity to pin 5 - Alternately switches with pin 1.      |
| 9                              | Nurse call relay - common connection          | ----                                                                                                 |

## Data transfer

### Upgrading software



**Upgrade of the Alaris® Syringe Pump (except the Alaris® TIVA Syringe Pump) software to V1.5.10 or greater is recommended at the next service, for Alaris® Syringe Pumps with software versions V1.5.9 and below. Perform upgrades by acquiring the software upgrade kits specified in spare parts listings. This will address a potential issue that may result in a condition where the running LED is flashing, the infusion status shows "INFUSING" but the volume infused display will not increment and no drug will be infused into the patient.**

**This potential issue may occur under the following remote circumstances :-**

- (i) A new syringe was recently fitted into the drive mechanism and
- (ii) An infusion is started, very quickly stopped and then restarted. (The pump must be stopped between 0.375 secs and 0.435 secs after starting - a window of 0.06 secs.)

**If the pump is subsequently stopped and restarted, the infusion will start normally.**

**When upgrading a pump from one software version to another where the first or middle digit changes, cold start will be required before and after software upgrade. Calibration will also be required after software upgrade and cold start.**

#### Tools required

- The Software Distribution Disk (See table below)
- Programming kit 1000SP00172 (Includes Programme Header and IrDA cable)
- IrDA port on PC or Comms Port
- RS232 cable 1000SP00336
- Ver. 3 Software Maintenance Utility (SMU) 1000CD00028

#### IrDA power-down test

To check PC is set up correctly for communication with Alaris® Syringe Pumps the Power Down Test needs to be performed on one Alaris® Syringe Pump only as follows:

1. Load the IrDA Power Down Test program on your PC.
2. Select GO on the PC software program.
3. Align the IrDA converter with the pump IrDA window (optimum distance is 5cm).
4. Connect to serial port.
5. Enter access code **166**.
6. Press Yes to continue Bootstrap.
7. Select IrDA interface.
8. Select a Baud rate of 115200.
9. The pump will then display Bootstrap in progress.
10. Press the  button to silence the alarm.
11. Select Transmit on PC. Check progress bar moves on PC and pump powers down.

#### Software Versions available

| Syringe Pump Model        | Software                 |                         | Enhanced Software       |                          | Guardrails® Safety Software                           |
|---------------------------|--------------------------|-------------------------|-------------------------|--------------------------|-------------------------------------------------------|
|                           | Mk1/Mk2                  | Mk3                     | Mk1/Mk2                 | Mk3                      | Mk3                                                   |
| Alaris® GS Syringe Pump   | 1000SP01221 (MP v1.5.10) | 1000SP01225 (MP v2.0.0) | 1000SP01270 (MP v1.9.4) | 1000SP01276 (MP v2.3.6)  |                                                       |
| Alaris® GH Syringe Pump   | 1000SP01221 (MP v1.5.10) | 1000SP01226 (MP v2.0.0) | 1000SP01270 (MP v1.9.4) | 1000SP01268 (MP v2.3.6)  | MP v3.1.4<br>(Installed by Cardinal Health Personnel) |
| Alaris® CC Syringe Pump   | 1000SP01221 (MP v1.5.10) | 1000SP01227 (MP v2.0.0) | 1000SP01270 (MP v1.9.4) | 1000SP01267 (MP v2.3.6)  | MP v3.1.4<br>(Installed by Cardinal Health Personnel) |
| Alaris® TIVA Syringe Pump | 1000SP01221 (MP v1.6.2)  | 1000SP01228 (MP v2.1.0) | 1000SP01270 (MP v1.9.4) | 1000SP01269 (MP v2.3.6)  |                                                       |
| Alaris® PK Syringe Pump   |                          |                         |                         | 1000SP01289 (MP v3.2.12) |                                                       |

Key: MP = Main Processor. Mk1/Mk2/Mk3 are the released versions of the Control PCB.

## Data transfer (continued)

### Soft bootstrap

1. Load the software program onto your PC. Start the 'MP Only' version of relevant pump software. Check the correct pump type is displayed.
2. Select GO.
3. Align the IrDA converter pump with the IrDA window (optimum distance is 5cm), or connect RS232 cable.
4. Connect to serial port.
5. Enter access code **166**.
6. Press Yes to continue Bootstrap.
7. Select IrDA interface or RS232 interface.
8. Select a Baud rate of 115200.
9. The pump will then display Bootstrap in progress.
10. Press the  button to silence the alarm.
11. Select Start on PC. Monitor progress of all selected channels
12. Power down pump.

### Hard bootstrap

1. Load the software program onto your PC. Start the relevant pump software (not the 'MP Only' version).
2. Disconnect the battery and separate the pump.
3. Fit the Programme header onto the control board.
4. Reconnect the battery. The pump will alarm, press the  button to silence.
5. Align the IrDA converter pump with the IrDA window (optimum distance is 5cm), or connect RS232 cable.
6. Connect to serial port.
7. Switch the Programme header to the correct position either RS232 or IrDA.
8. Switch on the Programme header.
9. Select GO on the PC software program.
10. Select Start on PC. Monitor progress of all selected channels
11. Power down pump.

### Cold start

It may be necessary to carry out a cold start if the pump has changed between certain software. Refer to documentation supplied with the software disk to see if cold start is required.

1. Enter access code **611**, then power down when prompted.
2. Perform a full calibration.



#### **Caution - Potential Erasure of Data:**

**Cold Start erases ALL information from the pump. This feature should only be used when changing between incompatible software versions. Full recalibration and reconfiguration will be required. Cardinal Health technicians should not re-instate drug information (this MUST be left to the customer).**



#### **Power Failure**

**Failures may occur when using laptops when communicating with Alaris® Syringe Pumps, due to power requirements. External power supply may be used in conjunction with IrDA or RS232 cable to compensate for lack of power from laptop.**

**Please Note IrDA data transfer can be affected by bright sunlight or fluorescent lighting.**

## Data transfer (continued)

### Teach Learn (Software versions V1.4.13 and above)

1. For the teacher pump only (not required for learn pumps), in General Options **251**, switch off RS232 before commencing Teach Learn procedure.
2. Turn the teacher pump on in normal operation.



**For multiple teach-learn operations, to avoid call-back alarm every 2 minutes, turn teacher pump on in access code mode.**

3. Enter the access code **167** on learn pump.
4. Align the two IrDA ports on the pumps (optimum distance is 5cm).
5. Press OK and then Yes to confirm.
6. A progress bar will travel across the learn pump.
7. When completed, select No to cancel retry.

#### Possible reasons for failure:

- RS232 is not switched off.
- If software versions are different, confirm teach learn procedure on learner pump to continue. Verify settings after completion of learn.
- The pump models are different. For example, an Alaris® GS Syringe Pump can only teach an Alaris® GS Syringe Pump.
- The line of sight between the IrDA windows was obstructed during data transfer.



**Check protocols are correct on learn pump after teach learn procedure, before returning pump to use.**

**After a Teach/learn procedure it is necessary to clear the previous patient setup in order to update the syringe information. On power-up after Teach/learn and when prompted with CLEAR SETUP, select YES.**

### Event Log download

A PC application known as the Event Log Download Utility (ELDU) (part number 1000SP00209) is available to download logs from Alaris® Syringe Pumps.

#### ELDU operation

1. Click on ELDU icon on PC.
2. Click Accept to agree with Restrictions of Use and continue.
3. Select Configure from drop-down menu.
4. Select Setup Pump and choose Alaris® as pump type.
5. Select Settings to select log to be downloaded.
6. Check communications are set as follows:
  - Required PC com port selected.
  - Set baud rate to 38400.
7. Click OK to confirm.
8. Align the IrDA converter pump with the IrDA window (optimum distance is 5cm), or connect RS232 cable.
9. Power up pump.
10. Click Download Log from main screen.
11. Press Close, when finished.
12. Select File from drop-down menu and save file. Log may be printed here as required.

## Data transfer (continued)

### Data Set Upload and Download (401 & 499)

#### Upload Data Set to an Alaris® Syringe Pump with Guardrails® Safety Software or an Alaris® PK Syringe Pump (401)

Using the Guardrails® Editor Transfer Tool or Alaris® PK Editor Software Transfer Tool allows a released Data Set to be uploaded to an Alaris® Syringe Pump.

#### Download Data Set from an Alaris® Syringe Pump with Guardrails® Safety Software or an Alaris® PK Syringe Pump (499)

Using the Verification Tool allows an uploaded Data Set in an Alaris® Syringe Pump to be downloaded to a PC for comparison and verification.

### Download CQI Event Log (402)

#### Download CQI Event Log from an Alaris® Syringe Pump with Guardrails® Safety Software (402)

Using the CQI Event Log Downloader allows the CQI Event Log to be downloaded from an Alaris® Syringe Pump to a PC for use with the Guardrails® CQI Reporter. The Guardrails® CQI Reporter is a program for querying and reporting on the collective event data allowing the user to analyse trends in medication administration and track medication errors.



#### Warning -

**At no time should the Guardrails® Safety Software or the Alaris® PK Editor Software be used to upload to or download from an Alaris® Syringe Pump while the pump is connected to a patient.**



**For more information relating to the Guardrails® Editor, the Alaris® PK Editor Software and the Guardrails® CQI Reporter refer to the relevant Directions For Use supplied with the software.**

## Information logs (376)

Use access code **376** to view the information logs (see Access Codes in chapter 2).

| Log           | View                                                                                                                                | Notes                                                         |
|---------------|-------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------|
| Service       | Displays the last 10 fault codes.                                                                                                   | Option to view the time and date at which they occur.         |
| Clear Service | Clears any information stored in the service log.                                                                                   | Will not be available if there is no data in the service log. |
| Event         | Displays the complete event log (maximum 1500 events except Pumps with Guardrails® Software enabled which have one year of events). | Option to view the time and date at which they occur.         |
| Key           | Displays the last 200 key presses and the time they occurred.                                                                       | Does not record while in Tech mode.                           |
| Use           | Displays the hours of use since reset and since last cold start.                                                                    | Press OK to clear hours since reset.                          |

Access code **376** provides the following additional service options:

- Service Date** Set the date when pump will display 'Service due' and any service message entered.
- Service Message** Enter message to be displayed on service date.
- Serial Number** Record the serial number of the pump.
- Unit Reference** Free-form text field for user reference only.
- Event Log** Access provided when standard power-up mode leads to errors such that the Event Log access from the Options ⓘ button cannot be accessed.
- PCB Identification Number** Allows Control PCB ID to be reviewed. (Pumps with Guardrails® Safety Software only)

### Recommended cleaning and inspection

To ensure this pump remains in good operating condition, it is important to keep it clean and carry out the routine procedures described below. All servicing should only be performed by a qualified service engineer.

- Thoroughly clean external surfaces of the pump, by wiping over with a lint-free cloth, lightly dampened with warm water and a standard disinfectant/detergent solution.

Do not use the following disinfectant types:

- NaDcc (such as PRESEPT)
- Hypochlorites (such as CHLORASOL)
- Aldehydes (such as CIDEX)
- Cationic Surfactants (such as Benzalkonium Chloride)
- Iodine (such as Betadine)

Recommended cleaners are:

| Brand     | Concentration |
|-----------|---------------|
| Hibiscrub | 20% (v/v)     |
| Virkon    | 1% (w/v)      |



***Before cleaning always switch OFF and disconnect from the AC power supply. Never allow fluid to enter the casing and avoid excess fluid build up on the pump.***  
***Do not use aggressive cleaning agents as these may damage the exterior surface of the pump.***  
***Do not steam autoclave, ethylene oxide sterilise or immerse this pump in any fluid.***

- Labels should be replaced as required if not flat, legible or fully adhered.
- Case components must be checked for damage and replaced if necessary.
- Check the pole clamp is not damaged and that it functions correctly.
- Inspect the AC power supply plug and cable for damage.



***Use an appropriate cleaning method that does not allow an excess of fluid to accumulate around the keypads. Aggressive cleaning can potentially create a fluid ingress path into the shelf keypad which can result in keypad failure.***

***In case of failure, usually resulting in a KY1 error code, the shelf keypad must be replaced. As a preventive measure, shelf keypads manufactured after week number 15, 2003 should be used since they offer more protection to excessive cleaning. The week number may be found on the keypad connection tail.***

***We recommend that all pumps within the following serial numbers -***

|                           |               |
|---------------------------|---------------|
| Alaris® GS Syringe Pump   | 08510 - 09976 |
| Alaris® GH Syringe Pump   | 16437 - 22286 |
| Alaris® CC Syringe Pump   | 03471 - 06632 |
| Alaris® TIVA Syringe Pump | 01310 - 02369 |

***(or pumps outside of this range which had their shelf keypad replaced between 2nd July 2002 and 30th April 2003) have their shelf keypad replaced at the next routine service. All other pumps have a shelf keypad that does not exhibit this potential risk.***

## Performance verification procedure

|                                               |                                                                                   |                                          |                          |                               |
|-----------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------|--------------------------|-------------------------------|
| <b>Model / Serial Number:</b>                 |                                                                                   | <b>Service Order / Inventory Number:</b> |                          |                               |
| <b>Hospital Name / Reference:</b>             |                                                                                   |                                          | <b>Software Version:</b> |                               |
| <b>INSPECTION</b>                             | <b>Physical inspection and clean</b>                                              |                                          |                          |                               |
| <b>UPDATES</b>                                | <b>Recommended when serviced</b>                                                  | UPDATE REF:                              | Fitted ✓                 | Not fitted / Not Applicable ✓ |
|                                               | For models GS, GH & CC upgrade software version to 1.5.10 or higher               | TSM <sup>CH3</sup>                       |                          |                               |
|                                               | Fit transmission buffer pad below 8001- 03469 & 8002 -06789                       | TSM <sup>CH6</sup>                       |                          |                               |
|                                               | Fit modified pole clamp arm as required                                           | TSM <sup>CH6</sup>                       |                          |                               |
|                                               | Fit support bracket for motor plate (1000SP00408)                                 | TSM <sup>CH6</sup>                       |                          |                               |
| <b>ERROR LOG</b> <sup>CH3</sup><br><b>376</b> | Check/set serial number, set service date (optional)                              |                                          |                          |                               |
| <b>SELF TEST</b> <sup>CH3</sup><br><b>123</b> | <b>Check all functions in self-test</b>                                           |                                          |                          |                               |
|                                               | Check date and time is correct (set as required (251) <sup>CH2</sup> )            |                                          |                          |                               |
| <b>INFUSING</b>                               | <b>Syringe size detection test</b>                                                |                                          |                          |                               |
|                                               | 12 mm spacer (11.5 to 12.5)                      32mm spacer (31.5 to 32.5)       |                                          |                          |                               |
| <b>VERIFICATION TESTS</b> <sup>CH3</sup>      | <b>Linear speed test*</b>                                                         |                                          | _____ mins _____ secs    |                               |
|                                               | <b>Occlusion test</b>                                                             |                                          | _____ KgF                |                               |
|                                               | <b>Line pressure readings (CC)</b>                                                |                                          | _____ mmHg               |                               |
|                                               |                                                                                   |                                          | _____ mmHg               |                               |
| <b>SETUP</b>                                  | <b>Set rate to zero (or lowest value possible), Clear Volume Infused and VTBI</b> |                                          |                          |                               |
| <b>ELECTRICAL SAFETY TESTS</b>                | <b>Class I Type CF</b> <i>Alternatively attach printed test results</i>           |                                          | _____ Ω                  |                               |
|                                               | Earth Resistance Test <= 0.2 Ω                                                    |                                          | _____ μA                 |                               |
|                                               | Earth Leakage Current <= 500 μA                                                   |                                          | _____ μA                 |                               |
|                                               | Enclosure Leakage Current <= 100 μA                                               |                                          |                          |                               |
| <b>Verification Performed By</b>              | _____                                                                             | _____                                    | _____                    |                               |
|                                               | <b>Sign</b>                                                                       | <b>Print</b>                             | <b>Date</b>              |                               |

<sup>CHX</sup> indicates the chapter number in the Technical Service Manual (TSM) - 1000SM00001.

E.G. <sup>CH3</sup> = Refer to TSM Chapter 3



**\* Latest issue of the plunger protector jig (0000JG00014 Issue 7) has been improved so that the needle of the dial gauge rests upon the plunger head (avoid resting the needle on the moulding flash line) of the pump. This improves the linear speed accuracy test results as any variation caused by the jig movement during the test are eliminated.**

## ***Troubleshooting***

### **In this chapter**

|                                 |           |
|---------------------------------|-----------|
| <b>Software fault codes</b>     | <b>31</b> |
| <b>Exception error handling</b> | <b>36</b> |
| <b>General fault diagnosis</b>  | <b>36</b> |

## Software fault codes



The following errors, MT1, DE1, PF1, PP1 and SC1 may be experienced if the self test operation or calibration operation has been accessed by quitting from the configuration menu. If these are displayed the pump should be power cycled and these operations entered directly.

| Code | Module              | Failure                                      | Action/Replace                                             |
|------|---------------------|----------------------------------------------|------------------------------------------------------------|
| AC1  | AC Alarm manager    | AC alarm manager failure                     | Control PCB                                                |
| AC2  |                     | AC VCO failure                               |                                                            |
| AD1  | ADC Converter       | Voltage reference/Power regulation           | Control PCB                                                |
| AM1  | Audio manager       | Audio status output driver                   | Control PCB                                                |
| AM2  |                     | VCO failure                                  |                                                            |
| AS1  | Audio Status        | Software execution                           | Speaker, wiring or Control PCB                             |
| AS2  |                     | Audio status monitoring input ADC            |                                                            |
| AS3  |                     | Speaker current test at power up             |                                                            |
| BT1  | Battery             | Battery gas gauge                            | Battery or Control PCB                                     |
| BT2  |                     | Battery cell voltage is low                  |                                                            |
| BT3  |                     | Battery cell voltage is high                 |                                                            |
| BT4  |                     | Battery discharging when connected to mains. |                                                            |
| CK1  | Clock               | Excessive timing drift                       | Control PCB                                                |
| DE1  | Drive Engage Detect | Drive engagement software module             | Control PCB                                                |
| DE2  |                     | Drive engagement opto self test              | Plunger drive flexi, Control PCB or Transmission PCB flexi |
| DE3  |                     | Emitter in wrong state                       |                                                            |
| DB1  | Drugs Manager       | Drug database file system                    | Control PCB                                                |
| DB2  |                     | Drug database file retrieval                 |                                                            |
| DB3  |                     | Drug database file storage                   |                                                            |
| DB4  |                     | CRC Error                                    |                                                            |
| DS1  | Dosing Manager      | Dosing retrieval failure                     | Control PCB                                                |
| DS2  |                     | Dosing storage failure                       |                                                            |
| DS3  |                     | Dosing data failure                          |                                                            |
| DS4  |                     | Dosing drug library failure                  |                                                            |
| DS5  |                     | Dosing patient data failure                  |                                                            |
| DS6  |                     | Dosing IDFS failure                          |                                                            |
| DS7  |                     | Dosing Data Set manager failure              |                                                            |
| DS8  |                     | Dosing Profile manager failure               |                                                            |
| DSM1 | Data Set Manager    | Data Set loading failure                     | Upload or check Data Set                                   |
| DSM2 |                     | Data Set integrity failure                   |                                                            |
| DSM3 |                     | Data Set incompatible format                 |                                                            |
| DSM4 |                     | Data Set CRC failure                         |                                                            |
| DSM5 |                     | Data Set read failure                        | Upload or check Data Set Control PCB                       |
| DSM6 |                     | Data Set element CRC failure                 |                                                            |
| DSM7 |                     | Data Set storage failure                     | Control PCB                                                |
| DSM8 |                     | Data Set storage retrieval failure           |                                                            |
| DSM9 |                     | Data Set data corruption                     |                                                            |

## Software fault codes (continued)

| Code | Module                 | Failure                             | Action/Replace                             |
|------|------------------------|-------------------------------------|--------------------------------------------|
| EV1  | Event Log              | Open log file at power up           | Control PCB                                |
| EV2  |                        | File storage software module        |                                            |
| EV3  |                        | Log read index                      |                                            |
| EV4  |                        | Log write index                     |                                            |
| EV5  |                        | Log data read                       |                                            |
| EV6  |                        | Log data write                      |                                            |
| EV7  |                        | Log data seek                       |                                            |
| EV8  |                        | Log repair failure                  |                                            |
| EV9  |                        | Log format failure                  |                                            |
| EV10 |                        | Log reporting failure               |                                            |
| EV11 |                        | Log extracting failure              |                                            |
| EV12 |                        | Log pack failure                    |                                            |
| EV13 |                        | Log unpack failure                  |                                            |
| FD1  | Fluid Delivery         | Fluid delivery software module      | Control PCB                                |
| FD2  |                        | Alarm manager software module       |                                            |
| FD3  |                        | Plunger drive software module       |                                            |
| FD4  |                        | Pressure monitor software module    |                                            |
| FD5  |                        | Syringe software module             |                                            |
| FD6  |                        | User configuration software module  |                                            |
| FD7  |                        | Fluid delivery setup data retrieval |                                            |
| FD8  |                        | Fluid delivery setup data storage   |                                            |
| FD9  |                        | Fluid delivery critical data        |                                            |
| FD11 |                        | VI cross check error                |                                            |
| FL1  | Fluid Log              | Open fluid log file at power up     | Control PCB                                |
| FL2  |                        | File storage software module        |                                            |
| FL3  |                        | Log read index                      |                                            |
| FL4  |                        | Log write index                     |                                            |
| FL5  |                        | Log data read                       |                                            |
| FL6  |                        | Log data write                      |                                            |
| FL7  |                        | Log data seek                       |                                            |
| FL8  |                        | Log repair failure                  |                                            |
| FL9  |                        | Log format failure                  |                                            |
| FL10 |                        | Log reporting failure               |                                            |
| FL11 |                        | Log extracting failure              |                                            |
| FL12 |                        | Log pack failure                    |                                            |
| FL13 |                        | Log unpack failure                  |                                            |
| GG1  | Gas Gauge              | Communications with gas gauge       | Battery or Control PCB                     |
| GR1  | Guardrails® Limit      | Guardrails® Limit failure           | Control PCB                                |
| IM1  | Identification Manager | Data corrupt                        | Control PCB                                |
| IM2  |                        | Serial number to set                | Set Serial Number in access code 376       |
| IM3  |                        | RTC to set                          | Set RTC in access code 251.<br>Control PCB |
| IM4  |                        | Storage failure                     | Control PCB                                |

**Software fault codes (continued)**

| Code | Module         | Failure                                                                  | Action/Replace                                                                                                            |
|------|----------------|--------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|
| KL1  | Key Log        | Open key log file at power up                                            | Control PCB                                                                                                               |
| KL2  |                | File storage software module                                             |                                                                                                                           |
| KL3  |                | Log read index                                                           |                                                                                                                           |
| KL4  |                | Log write index                                                          |                                                                                                                           |
| KL5  |                | Log data read                                                            |                                                                                                                           |
| KL6  |                | Log data write                                                           |                                                                                                                           |
| KL7  |                | Log data seek                                                            |                                                                                                                           |
| KL8  |                | Log repair failure                                                       |                                                                                                                           |
| KL9  |                | Log format failure                                                       |                                                                                                                           |
| KL10 |                | Log reporting failure                                                    |                                                                                                                           |
| KL11 |                | Log extracting failure                                                   |                                                                                                                           |
| KL12 |                | Log pack failure                                                         |                                                                                                                           |
| KL13 |                | Log unpack failure                                                       |                                                                                                                           |
| KY1  | Keypad         | Keypad key stuck down for 10mins                                         | Keypad or Control PCB                                                                                                     |
| KY2  |                | Line failure                                                             |                                                                                                                           |
| LC2  | LCD            | LCD control parameter storage                                            | Control PCB or Display PCB                                                                                                |
| LC3  |                | LCD display memory test                                                  |                                                                                                                           |
| ME1  | Motor Encoder  | Motor encoder module software                                            | Motor encoder or Control PCB                                                                                              |
| ME2  |                | Overrun detected                                                         |                                                                                                                           |
| ME3  |                | Motor encoder interrupt service software run when encoder is disabled    |                                                                                                                           |
| MT1  | Stepper Motor  | Motor module software                                                    | Check motor wire connections, opto flag not slipping and encoder gear is not in line with opto. Chassis PCB, Control PCB. |
| MT2  |                | Encoder has failed, preventing motor software continuing to run          |                                                                                                                           |
| MT3  |                | Safety processor has failed, preventing motor software continuing to run |                                                                                                                           |
| MT4  |                | Motor rotation in wrong direction                                        |                                                                                                                           |
| MT5  |                | Motor rotation speed has drifted                                         |                                                                                                                           |
| MT6  |                | Running at wrong rate                                                    |                                                                                                                           |
| MT7  |                | Motor rotation detected when it should be stopped                        |                                                                                                                           |
| MT8  |                | Motor not rotating when it should be                                     |                                                                                                                           |
| MT9  |                | Motor rotation inhibit control of Safety Processor                       |                                                                                                                           |
| MT10 |                | Motor rotation travel data has reached maximum value                     |                                                                                                                           |
| MT11 |                |                                                                          |                                                                                                                           |
| MT12 |                | Motor critical data                                                      |                                                                                                                           |
| NC1  | Nurse Call     | ADC preventing nurse call operation                                      | Control PCB or RS232/Nurse Call PCB                                                                                       |
| NC2  |                | Monitor signal is out of range                                           |                                                                                                                           |
| NC3  |                | Relay current monitor drive                                              |                                                                                                                           |
| NC4  |                | Safety processor failure                                                 |                                                                                                                           |
| PA1  | Patient Data   | Patient Data retrieval failure                                           | Control PCB                                                                                                               |
| PA2  |                | Patient Data storage failure                                             |                                                                                                                           |
| PA3  |                | Patient Data data error                                                  |                                                                                                                           |
| PA4  |                | Patient Data IDFS error                                                  |                                                                                                                           |
| PB1  | Plunger Button | Plunger button module software                                           | Plunger drive flexi, Control PCB or Transmission PCB flexi                                                                |
| PB2  |                | Plunger button opto self test                                            |                                                                                                                           |
| PB3  |                | Emitter in wrong state                                                   |                                                                                                                           |

## Software fault codes (continued)

| Code | Module                   | Failure                             | Action/Replace                                                                             |
|------|--------------------------|-------------------------------------|--------------------------------------------------------------------------------------------|
| PC1  | PIP Controller           | PIP Controller software module      | Control PCB                                                                                |
| PC2  |                          | PIP Controller critical data        |                                                                                            |
| PC3  |                          | Phase module critical data          |                                                                                            |
| PC4  |                          | Fluid delivery software module      |                                                                                            |
| PC5  |                          | PIP Controller setup data retrieval |                                                                                            |
| PC6  |                          | PIP Controller setup data storage   |                                                                                            |
| PD1  | Pressure Disc            | Pressure disc module software       | Pressure disc opto or Control PCB                                                          |
| PD2  |                          | Pressure disc opto self test        |                                                                                            |
| PD3  |                          | Emitter in wrong state              |                                                                                            |
| PF1  | Plunger Fitment          | Plunger fitment software module     | Optos,Cables, Control PCB or Plunger assy.                                                 |
| PF2  |                          | Gripper opto module software        |                                                                                            |
| PF3  |                          | Plunger button opto module software |                                                                                            |
| PFM1 | Profile Manager          | Profile storage failure             | Control PCB                                                                                |
| PFM2 |                          | Profile retrieval failure           |                                                                                            |
| PFM3 |                          | Profile startup failure             |                                                                                            |
| PFM4 |                          | Profile data corruption             |                                                                                            |
| PG1  | Plunger Grippers         | Plunger gripper module software     | Plunger gripper opto, Transmission PCB, or Control PCB                                     |
| PG2  |                          | Plunger gripper opto self test      |                                                                                            |
| PG3  |                          | Emitter in wrong state              |                                                                                            |
| PL1  | Plunger Drive            | Plunger drive software module       | Control PCB or Chassis PCB                                                                 |
| PL2  |                          | Alarm manager software module       |                                                                                            |
| PL3  |                          | Plunger drive travel deviation      | Calibrate linear travel. Plunger position linear potentiometer. Control PCB or Chassis PCB |
| PL4  |                          | Plunger position monitor software   | Control PCB or Chassis PCB                                                                 |
| PL5  |                          | Motor software module               |                                                                                            |
| PL6  |                          | Syringe software module             |                                                                                            |
| PL7  |                          | User configuration option software  |                                                                                            |
| PM1  | Pressure Measurement     | Pressure measurement software       | Control PCB                                                                                |
| PM2  |                          | Pressure sensor                     |                                                                                            |
| PM3  |                          | Syringe force                       |                                                                                            |
| PM4  |                          | User configuration options          |                                                                                            |
| PM5  |                          | Syringe                             |                                                                                            |
| PM6  |                          | Setup retrieval failure             |                                                                                            |
| PM7  |                          | Setup storage failure               |                                                                                            |
| PP1  | Plunger Position Monitor | Plunger position monitor software   | Control PCB                                                                                |
| PP2  |                          | ADC                                 |                                                                                            |
| PP3  |                          | Plunger position sensor             | Plunger position linear potentiometer.                                                     |
| PP4  |                          | Plunger position sensor calibration | Calibrate linear travel                                                                    |
| PP5  |                          | Plunger position cal data retrieval | Control PCB                                                                                |
| PP6  |                          | Plunger position cal data storage   |                                                                                            |
| PP7  |                          | Drive engagement software           |                                                                                            |

## Software fault codes (continued)

| Code | Module           | Failure                              | Action/Replace                             |
|------|------------------|--------------------------------------|--------------------------------------------|
| PS1  | Pressure Sensor  | Pressure sensor software module      | Control PCB                                |
| PS2  |                  | ADC                                  |                                            |
| PS3  |                  | Current reading invalid              | Pressure sensor or Control PCB             |
| PS4  |                  | Voltage (Normal) reading invalid     |                                            |
| PS5  |                  | Voltage (Test) reading invalid       |                                            |
| PS6  |                  | Pressure sensor calibration          | Calibrate pressure sensor                  |
| PS7  |                  | Pressure sensor cal data retrieval   | Control PCB                                |
| PS8  |                  | Pressure sensor cal data storage     |                                            |
| PS9  |                  | Pressure sensor amplifier gain       | Pressure sensor or Control PCB             |
| PS10 |                  | Pressure sensor shift failure        |                                            |
| SC1  | Syringe Clamp    | Syringe clamp software module        | Control PCB or syringe clamp potentiometer |
| SC2  |                  | ADC                                  |                                            |
| SC3  |                  | Syringe clamp sensor readings        |                                            |
| SC4  |                  | Syringe clamp calibration            | Calibrate syringe clamp                    |
| SC5  |                  | Syringe clamp cal data retrieval     | Control PCB                                |
| SC6  |                  | Syringe clamp cal data storage       |                                            |
| SD1  | SD Data          | Service data file retrieval          | Perform Cold start. Control PCB            |
| SD2  |                  | Service data file storage            |                                            |
| SD3  |                  | Service data file contents           |                                            |
| SF1  | Syringe Force    | Syringe force sensor software        | Control PCB                                |
| SF2  |                  | ADC                                  |                                            |
| SF3  |                  | Syringe force sensor current reading | Motor plate or Control PCB                 |
| SF4  |                  | Syringe force sensor normal reading  |                                            |
| SF5  |                  | Voltage (test) reading               |                                            |
| SF6  |                  | Syringe force sensor calibration     | Calibrate syringe force                    |
| SF7  |                  | Syringe force cal data retrieval     | Control PCB                                |
| SF8  |                  | Syringe force cal data storage       |                                            |
| SF9  |                  | Syringe force sensor amplifier       | Motor plate or Control PCB                 |
| SF10 |                  | Syringe force output                 | Calibrate syringe force. Motor plate.      |
| SP1  | Safety Processor | Safety processor program memory      | Control PCB                                |
| SP2  |                  | Safety processor ID version          |                                            |
| SP3  |                  | Safety processor                     |                                            |
| SV1  | Service Log      | Service Log software module          | Control PCB                                |
| SV2  |                  | File storage software module         |                                            |
| SV3  |                  | Log read index                       |                                            |
| SV4  |                  | Log write index                      |                                            |
| SV5  |                  | Log data read                        |                                            |
| SV6  |                  | Log data write                       |                                            |
| SV7  |                  | Log data seek                        |                                            |
| SV8  |                  | Log repair failure                   |                                            |
| SV9  |                  | Log format failure                   |                                            |
| SV10 |                  | Log reporting failure                |                                            |
| SV11 |                  | Log extracting failure               |                                            |
| SV12 |                  | Log pack failure                     |                                            |
| SV13 |                  | Log unpack failure                   |                                            |

## Software fault codes (continued)

| Code   | Module               | Failure                            | Action/Replace                                                              |
|--------|----------------------|------------------------------------|-----------------------------------------------------------------------------|
| SY1    | Syringe Manager      | Syringe manager software module    | Control PCB                                                                 |
| SY2    |                      | Syringe clamp                      |                                                                             |
| SY4    |                      | Plunger fitment                    |                                                                             |
| SY5    |                      | Plunger button stuck in            | Check plunger button. Control PCB.                                          |
| SY6    |                      | Syringe data table retrieval       | Control PCB                                                                 |
| SY7    |                      | Syringe data table storage         |                                                                             |
| TC0    | PK / TCI             | No fault found                     | Contact your local Cardinal Health, Alaris® Products service representative |
| TC1    |                      | TCI data corruption                |                                                                             |
| TC2    |                      | TCI not configured                 |                                                                             |
| TC3    |                      | TCI incorrect state                |                                                                             |
| TC4-12 |                      | PK Model error 4-12                |                                                                             |
| TC13   |                      | PK Model parameters error          |                                                                             |
| TC14   |                      | PK Model init error                |                                                                             |
| UC1    | User Config. options | User config. option file retrieval | Control PCB                                                                 |
| UC2    |                      | User config. option file storage   |                                                                             |
| UC3    |                      | User config. option data retrieval |                                                                             |
| UC4    |                      | User configuration model           | Perform Cold start                                                          |
| UT1    | User Timer           | User timer master clock            | Control PCB                                                                 |
| UT2    |                      | User timer software request        |                                                                             |

## Exception error handling

Exception errors include Assertion Errors and Enum Failure Errors and are used to trap logical errors in the software execution.

The pump will display the error type, the title of the software module in which the error occurred and the line number. The user should make a note of these for use in diagnosis. This information is stored in the service log (access code **376**).

After an error, the pump will not store information when powered down. When the pump is switched on again, the user should always confirm clear setup (this is done automatically with software version V1.5.9 and later).

## General fault diagnosis


|               |                                 | Parts to Check/Test |           |                  |           |             |           |             |         |            |       |
|---------------|---------------------------------|---------------------|-----------|------------------|-----------|-------------|-----------|-------------|---------|------------|-------|
|               |                                 | Front Case          | Rear Case | Labels & Keypads | Mechanism | Control PCB | Power PCB | Display PCB | Battery | Mains Lead | Fuses |
| General Fault | Dropped or damaged              | ✓                   | ✓         |                  | ✓         | ✓           | ✓         | ✓           |         |            |       |
|               | Exposed to fluids               | ✓                   | ✓         | ✓                | ✓         | ✓           | ✓         | ✓           | ✓       | ✓          | ✓     |
|               | No battery power                |                     |           | ✓                |           | ✓           |           |             | ✓       |            |       |
|               | No AC mains power               |                     |           | ✓                |           | ✓           | ✓         |             |         | ✓          | ✓     |
|               | Delivery rates out of tolerance | ✓                   |           |                  | ✓         | ✓           |           |             |         |            |       |

## ***Circuit Descriptions***

### **In this chapter**

|                                               |           |
|-----------------------------------------------|-----------|
| <b>Functional module block diagram</b>        | <b>38</b> |
| <b>Module overview functional description</b> | <b>39</b> |

## Functional module block diagram



## Module overview functional description

The Alaris® Syringe Pumps are designed to be serviced generally to major assembly level. The PCBs are designed as non-serviceable items and as such, can only be replaced as complete parts.

The major assemblies are:

- Control PCB
- Power Supply Unit PCB
- Display PCB
- Battery
- Transmission
- Transducer (Model CC)

Cardinal Health will make available, on request, circuit diagrams which will assist appropriately qualified technical personnel to repair those parts of the device which are designated by the manufacturer as repairable.

### Control PCB

Contains the main processor module, which provides the control functions for almost all aspects of the pump. It drives and monitors all other modules using the program code stored in the flash eprom. The main processor runs the main application program. The main processor directly interfaces to:

- Safety Processor
- Keypads
- Display
- Real Time Clock
- Communications Switch to IrDA and RS232 (optional) interfaces
- Nurse Call Output
- Indicator LEDs
- Audible Alarm
- Motor – controller and dual channel coil driver DAC
- System Sensors (including: syringe clamp, plunger position, drive engagement opto, plunger button opto, syringe force sensor, line pressure sensor, pressure disc opto, motor encoder, AC power).
- Backlight
- Power supply
- Battery gas gauge

The function of the Safety Processor Module is to ensure the correct operation of the Main Processor by constantly exchanging data. If an error is detected, the module can independently disable the stepper motor that drives the transmission. Additionally, it can create both audible and visible alarms using its dedicated piezoelectric buzzer, alarm LED and, if fitted, the Nurse Call Interface.

The Safety Processor controls (independent of main processor):

- Audio sounder
- Visual indicator LED
- Control signal to inhibit motor drive
- Power supply hold up control

### Pressure Transducer (Model CC)

Monitors the line pressure when the pressure disc is inserted and flags the presence of a pressure disc. The Control PCB checks the transducer for presence of a pressure disc and the line pressure when a disc is present.

### Power Supply Unit PCB

The Power Supply Unit (PSU) is a single output switched-mode power supply rated at 15 VDC/1.6 A continuous duty. The PSU has a wide input voltage range of 85 to 264VAC, 47-63 Hz single phase.

### Display PCB

The Pump uses a Cold Cathode Fluorescent Lamp (CCFL) as a backlight for the negative mode LCD display. The CCFL Backlight Supply Module generates the high voltages required to drive the lamp and facilitates software based brightness control.

### Battery

The Battery Pack Module provides system power in the absence of a mains supply.

The Battery Pack Module consists of six 1.2V 2.7Ah NiMH battery cells connected in series, a thermal fuse, thermal circuit breaker and Gas Gauge Module sealed in a heat shrink sleeving.

The Gas Gauge Module is permanently connected across the battery terminals so that it can monitor terminal voltage, charge / discharge current and the battery pack temperature.

Through charge monitoring information, from the Gas Gauge Module, the Control PCB Main Processor Module determines the battery charge level and hence 'Battery Low' and 'Battery Empty' conditions.

Battery capacity will reduce over time.

Annual battery calibration is recommended to update Gas Gauge Module with 'up to date' battery capacity information.

If the battery pack fails to achieve the calibration limits, it is recommended that the battery pack is replaced and calibration performed.

### Transmission

The Transmission Interface Module provides the Pump with the capability of monitoring a number of critical parameters associated with the transmission operation. The device can detect failures or incorrect operation of the transmission and prevent incorrect drug dosage administration.

The electronics for this module occupy three separate PCBs that are located in various areas of the transmission as follows:

**Plunger PCB** – Located inside the Plunger Holder Assembly. Contains the electronics that detect correct syringe plunger location and gripper motion.

**Carriage PCB** – Located on the Transmission Carriage. Contains electronics that detect correct engagement of the Half Nut with the Leadscrew.

**Chassis PCB** – Located on the Chassis Extrusion. This facilitates measurement of motor speed, syringe drive force and linear plunger position. Additionally, the IrDA compliant infrared transceiver is positioned on the reverse side of this PCB. The stepper motor which drives the mechanical transmission of the pump is controlled by the Motor Drive module on the Control PCB.

## ***Spare Parts Replacement Procedures***

### **In this chapter**

|                                                |           |
|------------------------------------------------|-----------|
| <b>Access to pump</b>                          | <b>41</b> |
| <b>Rear case and subassemblies</b>             | <b>43</b> |
| <b>Front case and subassemblies</b>            | <b>47</b> |
| <b>Pressure transducer assembly (Model CC)</b> | <b>57</b> |
| <b>Keypads and labels</b>                      | <b>58</b> |

## Spare Parts Replacement Procedures



Ensure the pump is disconnected from the AC power supply and switched off before attempting to service. The pump contains static-sensitive components and therefore strict ESD precautions should be observed at all times.

Always protect the plunger holder and syringe clamp when the pump is upside down. For regular servicing, the use of the case support cradle Part No. 0000JG00047 is recommended.

Batteries should be disposed of as outlined by the local country regulations. Do not send batteries back to the manufacturer.

For fastener torque settings, please refer to Appendix D Fitting & Replacement guidelines.


Only use Cardinal Health recommended spare parts.

Following all spare part replacement and repair activities, testing must be performed in accordance with the Performance Verification Procedure (PVP), see Chapter 3, 'Routine Maintenance'.

### Access to pump

#### Replacement Procedure

1. Remove the two case screws in battery cover, remove cover and battery.
2. Remove the six case screws.
3. Model CC only: Insert a flat-blade screwdriver into the blanking plug of transducer, prise plug away from transducer and remove securing screw.
4. Carefully separate case halves and disconnect cables.
5. Where necessary, remove the foot rivets with a flat-blade screwdriver and remove the feet from the case. Refer to additional information on the following page concerning rivet orientation.
6. Reassemble in reverse order.



#### Spare Parts

| Item | Description                                            | Part Number |
|------|--------------------------------------------------------|-------------|
| A    | ASENA SP, Assy, Battery                                | 1000SP01122 |
| B    | ASENA SP, Battery Cover/Handle                         | 1000SP01121 |
| C    | ASENA CC, Assy, Plug Blanking Transducer               | 1000ME01317 |
| D    | ASENA GS, Kit, Front Case                              | 1000SP00478 |
| D    | ASENA GH, Kit, Front Case                              | 1000SP00479 |
| D    | ASENA TIVA, Kit, Front Case                            | 1000SP00480 |
| D    | ASENA CC, Kit, Front Case                              | 1000SP01153 |
| D    | ASENA PK, Kit, Front Case                              | 1000SP01204 |
| E    | ASENA GS/GH/TIVA, Kit, Rear Case                       | 1000SP01115 |
| E    | ASENA CC, Kit, Rear case                               | 1000SP01154 |
| F    | ASENA CC, Kit, Pressure Transducer                     | 1000SP01155 |
| G    | ASENA SP, Case Sealing Cord (1m) (internal, not shown) | 1000ME00311 |
| H    | ASENA SP, Kit, Fixings (screws, washers etc)           | 1000SP00466 |
| I    | ASENA SP, Kit, Spare adhesive foot rivet replacement   | 1000SP00593 |
| I    | ASENA SP, Kit, Spare adhesive foot                     | 1000SP00595 |

### **Access to pump** *(continued)*

The Pump has two types of rivet feet, one type is white and the other is black. The rivet feet can be removed and replaced without opening the pump, and may be replaced with adhesive backed feet as described below.

There are also two types of cases, one that had rivet feet and therefore holes in the cases for them and the other is without the holes but only a recess for the adhesive backed foot to fit into.

For cases that had rivet feet and therefore holes, kit 1000SP00593 is required. For cases without the holes, kit 1000SP00595 is required.

Kit 1000SP00593 contains a fitting instruction, 242 adhesive backed feet, 242 foot bonding pads and one 20g tube of Loctite 454 gel adhesive. This kit is intended as a replacement for the rivet style foot and has enough feet for 48 pumps.


Kit 1000SP00595 contains 11 adhesive backed feet only. This kit is intended for cases that do not have the rivet feet and therefore no holes in the cases. Only a recess is present for the adhesive backed foot to be placed into.

## Rear case and subassemblies

## Power Supply Unit &amp; Speaker

## Replacement Procedure

1. Disconnect the PSU cable.
2. Remove the three PSU screws.
3. Remove earth wire screw and washer.
4. Remove PSU and insulator.
5. With a pair of soft-faced pliers, carefully compress the catch holding the internal speaker and pull the speaker up and out.
6. Reassemble in reverse order.



## Spare Parts


| Item | Description                                  | Part Number |
|------|----------------------------------------------|-------------|
| A    | ASENA SP, Kit, Speaker                       | 1000SP01130 |
| B    | ASENA SP, Assy, Power Supply Unit (PSU)      | 8000EL00063 |
| C    | ASENA SP, Assy, PSU Insulator                | 1000ME01306 |
| D    | ASENA SP, Kit, Fixings (screws, washers etc) | 1000SP00466 |

## Rear case and subassemblies *(continued)*

### Rear case – mains inlet, PE stud and magnet

#### Replacement Procedure

1. Remove two nuts to remove PE stud.
2. Remove the two screws on Mains inlet.
3. Remove mains inlet and retainer.
4. Remove magnet by lifting one end.
5. Reassemble in reverse order.



#### Spare Parts

| Item | Description                                          | Part Number |
|------|------------------------------------------------------|-------------|
| A    | ASENA SP/GW, Kit, PE Stud                            | 1000SP00467 |
| B    | ASENA SP, Kit, Mains Inlet                           | 1000SP01124 |
| C    | Magnet IR Detect                                     | 1000ME01303 |
| D    | ASENA SP, Fuse, T-1.25A Slow Blow, Mains (not shown) | 1000EL00222 |
| E    | ASENA SP, Kit, Fixings (screws, washers etc)         | 1000SP00466 |



## Rear case and subassemblies (continued)

## Pole clamp and RS232


## Replacement Procedure

1. Remove three pole clamp screws.
2. Remove two tube restraint screws.
3. Remove two nuts & washers from RS232 socket screws.
4. Reassemble in reverse order.

(E) RS232 Socket Connector



**The Pole Clamp Arm material has been changed to a stronger material to prevent the arm from bending when tightened.**  
**The Pole Clamp Arm spares kit replaces parts of the Pole Clamp assembly to address bent or slipping Pole Clamps.**  
**Note: There is no requirement to remove the V Clamp.**



## Spare Parts


| Item | Description                                  | Part Number |
|------|----------------------------------------------|-------------|
| A    | ASENA SP, Assy, Pole clamp                   | 1000SP00115 |
| B    | ASENA SP, Tube restraint blank               | 1000ME01213 |
| C    | ASENA SP, Tube restraint RS232               | 1000ME01214 |
| D    | SPARE KIT POLE CLAMP ARM                     | 1000SP00589 |
| E    | ASENA SP/GW, Kit, RS232 connector            | 1000SP00468 |
| F    | ASENA SP, Kit, Fixings (screws, washers etc) | 1000SP00466 |
| G    | POLE CLAMP SNAKE EYE DRIVER (NOT SHOWN)      | 1000ME01466 |
| H    | Spring Pole Clamp                            | 0000ME00379 |
| I    | Bearing ball 5mm SS AISI316 Grade 100        | 0000ME00369 |
| J    | Clamp Pole Pivot (shoulder screw)            | 1000ME01244 |

## Rear case and subassemblies (continued)

### Rail cam

#### Replacement Procedure

1. Remove screw from lever release.
2. Remove screw from lever rail cam.
3. Remove locking washer from spring.
4. Reassemble in reverse order.



### Spare Parts


| Item | Description             | Part Number |
|------|-------------------------|-------------|
| A    | ASENA SP, Kit, Rail Cam | 1000SP01114 |

## Front case and subassemblies

### Front case – Control PCB and RS232 (if option fitted)

#### Replacement Procedure

1. Disconnect cable from RS232 PCB.
2. Remove the four retaining screws and washers from RS232 PCB.
3. Remove the three retaining screws and disconnect all flexi and cable connections. Push the motor towards the transmission to ease removal of Control PCB.
4. When fitting Control PCB ensure all flexi and cables are routed clear of PCB.
5. Connect all flexi and cable connections - secure with the three retaining screws.
6. Reassemble RS232 PCB in reverse order.



#### Spare Parts


| Item | Description                                  | Part Number |
|------|----------------------------------------------|-------------|
| A    | ASENA GS, Control PCB (Software V2.2.1)      | 1000SP00494 |
| A    | ASENA GH, Control PCB (Software V2.2.1)      | 1000SP00495 |
| A    | ASENA CC, Control PCB (Software V2.2.1)      | 1000SP00496 |
| A    | ASENA TIVA, Control PCB (Software V2.2.1)    | 1000SP00497 |
| A    | ASENA PK, Control PCB (Software V3.2.5)      | 1000SP00637 |
| B    | ASENA SP, Kit, RS232 (PCB & Fixings)         | 1000SP01160 |
| C    | ASENA SP, Kit, Fixings (screws, washers etc) | 1000SP00466 |

## Front case and subassemblies (continued)

### Display PCB

#### Replacement Procedure

1. Remove the flexible circuit ferrite.
2. Remove the four display fixing screws and two display mounting brackets.
3. Reassemble in reverse order. For the Model GS fit the adaptor bracket to the display.
4. Secure the shelf keypad flexi to the display using double-sided adhesive pad.



#### Spare Parts


| Item | Description                                  | Part Number                                       |
|------|----------------------------------------------|---------------------------------------------------|
| A    | ASENA GS, Display PCB                        | 1000SP01118                                       |
| B    | ASENA GH/CC/TIVA, Display PCB                | 1000SP01119                                       |
| C    | ASENA GS, Assy, Window Adaptor Display       | 1000SP00122                                       |
| D    | ASENA SP, Assy, Bracket Display Mounting     | 1000ME00261                                       |
| E    | Pad Self-adhesive Double-sided 12x12mm       | 0000ME00423                                       |
| F    | ASENA GS, Assy, Display Insulator            | (not shown on back of Display PCB)<br>1000SP00187 |
| G    | ASENA GH/CC/TIVA, Assy, Display Insulator    | 1000SP00188                                       |
| H    | ASENA SP, Assy, Gasket Display               | 1000ME01301                                       |
| I    | ASENA GH/CC/TIVA, Kit, Display Protection    | 1000SP01137                                       |
| J    | ASENA GS, Kit, Backlight Shield              | (not shown on back of Display PCB)<br>1000SP00273 |
| K    | ASENA SP, Kit, Fixings (screws, washers etc) | 1000SP00466                                       |

## Front case and subassemblies (continued)


## Front case – Chassis PCB and Plunger assembly

## Replacement Procedure

1. Remove the two Chassis PCB screws. Disconnect all cables.
2. Extend the plunger out to its full extent and fully loosen the two plunger retaining screws in the carriage.
3. Carefully remove the plunger flexi from the carriage PCB and straighten. While applying controlled force to the plunger, extract it from the carriage and withdraw it completely.
4. Reassemble in reverse order.



**Check Buffer Pad fitted if serial numbers are within either of the ranges 8001-03468 & below or 8002-06788 & below. If not fitted, clean the surface of the carriage face nearest the plunger drive tube and fit Buffer Pad in the position shown (sloping edge to match carriage profile, see diagram).**



## Spare Parts


| Item | Description                                  | Part Number |
|------|----------------------------------------------|-------------|
| A    | ASENA SP, Kit, Chassis PCB                   | 1000SP00189 |
| B    | ASENA SP, Kit, Plunger assembly (not shown)  | 1000SP01113 |
| C    | ASENA SP, Kit, Carriage buffer               | 1000SP00230 |
| D    | ASENA SP, Kit, Fixings (screws, washers etc) | 1000SP00466 |

## Front case and subassemblies (continued)

### Chassis assembly and Pressure Transducer (Model CC only)

#### Replacement Procedure


1. Carefully peel away the label on the bottom of the front case to gain access to the chassis screw. Remove this screw.
2. Remove the two screws securing the syringe flange clamp.
3. Carefully withdraw the chassis.
4. Remove the two screws from the pressure transducer assembly and carefully withdraw (Model CC only).
5. Reassemble in reverse order.



**Syringe Flange Clamp has been enhanced in response to market feedback indicating under certain conditions, false error alarms may occur on the pump close to the End of Infusion (EOI) particularly when syringes of small sizes are used (5,10ml etc.) if incorrectly fitted to the pump.**

**Fit Enhanced 1 piece Syringe Flange Clamp if serial numbers are within either of the ranges 8001-02315 & below or 8002-04311 & below.**

**Alternatively the latest 2 piece Syringe Flange Clamp may be fitted if the pump is required to have an EOI point below 5%. Pump must have software versions v1.8.1 or higher if fitting the 2 piece Syringe Flange Clamp.**



#### Spare Parts


| Item | Description                                    | Part Number |
|------|------------------------------------------------|-------------|
| A    | ASENA SP, Kit, Chassis Assembly                | 1000SP01112 |
| B    | ASENA SP, Kit, Syringe Flange Clamps (1 piece) | 1000SP00577 |
| B    | ASENA SP, Kit, Syringe Flange Clamps (2 piece) | 1000SP00570 |
| C    | ASENA CC, Kit, Pressure Transducer             | 1000SP01155 |
| D    | ASENA SP, LBL, Label Chassis Screw Cover       | 1000LB00431 |
| E    | ASENA SP, Kit, Fixings (screws, washers etc)   | 1000SP00466 |

## Front case and subassemblies (continued)

### Front case – Syringe Sizing assembly


#### Replacement Procedure

1. Remove the syringe sizing retainer screws, case brace and retainer.
2. Remove the shim and discard, then lever the syringe sizing mechanism from the housing and withdraw the potentiometer.
3. Pull the syringe clamp back to its full extent. Carefully remove the bung, screw and washer. Pull hard on the syringe clamp to remove.
4. Carefully lever the syringe sizing mechanism from the housing and pull through the case. Remove the shaft bearing and the v-ring seal.
5. Secure the assembly loading jig to the syringe sizing mechanism. Fit the shaft bearing and v-ring seal onto the end of the jig.
6. Lay the assembly on one side, potentiometer to the left, wires exiting upwards. The injection 'pip' feature on the pre-moulded shaft should be visible.



7. Fit the seal protector into the upper case and load the syringe sizing mechanism. Compress the v-seal against the protector.
8. Withdraw the protector and push the syringe sizing mechanism through the hole in the front case until the flat sides locate in the case and the potentiometer aligns with the case recess. Ensure the moulding pip is located on the side.
9. Slide the shim component down the side wall of the syringe potentiometer recess. Bend the shim 'outward'.
10. Fit the syringe sizing retainer so that the shim is visible protruding from the retainer. Fit the case brace. Secure with two screws.
11. Remove the assembly loading jig.
12. The syringe shaft flats to be moved into the open position.
13. Fit the syringe clamp over the shaft, fit the screw, washer and bung into the shaft end.

## Front case and subassemblies (continued)



**Check for presence of shim. If shim is not fitted, or the cam has sharp edges, fit the Cam Kit on reassembly. If the cam has sharp edges or the shim is folded incorrectly, these may cause excessive wear of areas around the front case. Mechanical movement and small changes in the syringe diameter can result in a syringe detect failure, which may occur if shim is not present or case is worn.**



**If pump has serial numbers within either of the ranges 8001-02574 & below or 8002-04778 & below fit enhanced syringe clamp 1000SP01123 on reassembly. Old syringe clamps can be recognised as clear plastic; new syringe clamps are solid blue plastic. Old syringe clamps may crack and fail after being subjected to Isopropyl alcohol used in the cleaning process.**



**Replace the foot last if replacing the case. If fitting a new syringe sizing mechanism, apply generous amounts of grease (CASTROL LMX) to the slot of the mechanism and lightly grease the v-seal.**

### Spare Parts


| Item | Description                                                | Part Number |
|------|------------------------------------------------------------|-------------|
| A    | ASENA SP, Kit, CAM Kit                                     | 1000SP00170 |
| B    | ASENA SP, Kit, Syringe Sizer (see previous page)           | 1000SP01116 |
| C    | ASENA SP, Kit, Syringe & Flange Clamps (see previous page) | 1000SP01123 |
| D    | ASENA Kit Grease LMX                                       | 1000SP01144 |
| E    | ASENA SP, Assy, Syringe Size Potentiometer                 | 1000SP01125 |
| F    | ASENA SP, Kit, Fixings (screws, washers etc)               | 1000SP00466 |
| G    | ASENA SP, Test, Syringe Clamp Jig                          | 1000SP00481 |

## Front case and subassemblies (continued)

### Chassis assembly breakdown

#### Replacement Procedure

- 1 Remove the pulley nut, washers and withdraw the pulley and toothed belt.
- 2 Remove three screws securing stepper motor.
- 3 Remove the leadscrew by driving out the roll pin using a suitable punch.
- 4 Remove three screws securing motor plate.
- 5 Remove three screws securing bearing block.
- 6 Refit plunger into carriage, declutch plunger and withdraw plunger and carriage. Hold linear potentiometer actuator and spring on the side of the carriage.
- 7 Remove one screw holding carriage PCB.
- 8 Remove the linear travel potentiometer.
- 9 Fit new linear potentiometer to centre area of chassis and flush to rear of chassis slot and flush to motor-plate end.
- 10 Reassemble in reverse order



## Front case and subassemblies (continued)

### Chassis assembly breakdown (continued)




Figure 1 - Current Motor Plate




Figure 2 - Motor Plate Strain Beam Support



**Check Motor Plate serial code, if the code begins with "PH"; e.g. PHDA 123, then this is the current version (see Figure 1) and does not require a Motor Plate Strain Beam Support. All previous version Motor Plates require a Motor Plate Strain Beam Support fitted in position shown (see Figure 2) if not fitted.**

### Spare Parts


| Item | Description                                  | Part Number |
|------|----------------------------------------------|-------------|
| A    | ASENA SP, Kit, Chassis Assembly              | 1000SP01112 |
| B    | ASENA SP, Kit, Stepper Motor                 | 1000SP01109 |
| C    | ASENA SP, Kit, Motor Plate                   | 1000SP01110 |
| D    | ASENA SP, Assy, Linear Potentiometer         | 1000ME01451 |
| E    | ASENA SP, Carriage PCB                       | 8000EL00022 |
| F    | ASENA SP, Kit, Chassis Enhancement           | 1000SP01136 |
| G    | ASENA SP, Kit, Fixings (screws, washers etc) | 1000SP00466 |
| H    | Spare Carriage Asena                         | 1000SP01107 |
| I    | SP Kit Motor Plate Strain Beam Support       | 1000SP00408 |

## Front case and subassemblies (continued)

### Plunger assembly breakdown

#### Replacement Procedure

1. Remove three screws holding plunger backplate.
2. Remove components as required. Follow diagrams step ii to step ix.
3. Reassemble in reverse order.




#### Spare Parts

| Item | Description                                                               | Part Number |
|------|---------------------------------------------------------------------------|-------------|
| A    | ASENA SP, Kit, Plunger Assembly<br>(Complete Plunger Assembly, not shown) | 1000SP01113 |
| B    | ASENA SP, Plunger Detect PCB                                              | 8000EL00019 |
| C    | ASENA SP, Gripper Bottom                                                  | 1000ME01218 |
| D    | ASENA SP, Gripper Top                                                     | 1000ME01219 |
| E    | ASENA SP, Kit, Fixings (screws, washers etc)                              | 1000SP00466 |
| F    | GEAR DECLUTCH                                                             | 1000ME01198 |

Front case and subassemblies (continued)

Plunger assembly breakdown (continued)




## Pressure Transducer Assembly (Model CC only)

The following instructions detail the fitting of the Pressure Transducer Assembly.

### Replacement Procedure


1. Fit the mylar gasket.
2. Align the hole in the gasket with the hole in the centre disc holder and ensure label is square to the Centre Disc Holder.
3. Use a clean wipe and apply pressure to the mylar gasket.
4. Crease the mylar gasket along the ledge of the centre disc holder and ensure it is well adhered along the front face of the step edge.
5. Load the spring onto the disc-detect flag shaft.
6. Locate spring arms to spring retainer on the flag and to the recess in the disc holder top.




13. Slide the flexible circuit ferrite into place and ensure the flexible circuit is formed at 90° to the PCB.
14. Lower the Base Disc Holder onto the Centre Disc Holder.
15. Secure the base with the screws.
16. Use a lint-free cloth and approved cleaner to wipe the surface clean.
17. Start the gasket at the pointed end. Work around the perimeter and minimise the gap at the join (if using cord - N/A if using single-piece gasket).

### Spare Parts

| Item | Description                                  | Part Number |
|------|----------------------------------------------|-------------|
| A    | ASENA CC, Kit, Pressure Transducer           | 1000SP01155 |
| B    | ASENA CC, Assy, Top Disc Holder              | 1000ME00450 |
| C    | ASENA CC, Assy, Flag Disc Holder             | 1000ME01318 |
| D    | ASENA CC, Assy, Gasket Mylar                 | 1000ME01319 |
| E    | ASENA SP, Kit, Fixings (screws, washers etc) | 1000SP00466 |



7. Rotate and install the plastic flag.
8. Fit the sealing cord starting at the break bar.
9. Load the disc holder centre onto the disc holder top.
10. Secure disc holder centre to disc holder top using screw.
11. Load pressure transducer into assembled housing.
12. Once the PCB is located, apply pressure over the sensor area of the PCB to ensure good location.




## Keypads and labels

### Replacement Procedure

1. Discard any keypads removed as they cannot be reused.
2. Ensure all residual adhesive is removed from bonding surfaces.
3. Fit replacement keypads after removing backing paper from underside. Handle replacement keypads carefully to avoid damage.
4. Apply finger pressure to keypads working from one end, to drive the air out of the adhesive/case interface.
5. Remove label(s) from case as required.
6. Clean case where replacement label(s) to be fitted.
7. Fit replacement label(s) taken from label sheet as required.
8. Ensure all keypad membrane flexi tails are routed and secured with double sided adhesive pads (0000ME00423).



**To ensure an effective case fluid seal, ensure the top 5mm edge of the shelf keypad is given careful attention as described in step 4 above.**




### Spare Parts

| Item | Description                             | Part Number |
|------|-----------------------------------------|-------------|
| A    | ASENA GH/CC, Key, Switch Overlay Keypad | 1000SP01126 |
| B    | ASENA GS, Key, Switch Overlay Keypad    | 1000SP01127 |
| C    | ASENA TIVA, Key, Switch Overlay Keypad  | 1000SP00403 |
| D    | ASENA PK, Key, Switch Overlay Keypad    | 1000SP01203 |


## Keypads and labels (continued)

### Note:

Item **A** additionally provides a clear window for the combined serial number & status label part nr. 1000LB00102.



**Note:** Item **B** also available separately. Part No. 1000LB00431



Picture shows Alaris® GS Syringe Pump Label Set. Refer to the following model types for correct label set.


### Description

LABEL SET CC SYRINGE PUMP  
 LABEL SET GH SYRINGE PUMP  
 LABEL SET GS SYRINGE PUMP  
 LABEL SET TIVA SYRINGE PUMP  
 LABEL SET PK SYRINGE PUMP


### Part Number

1000LB01408  
 1000LB01409  
 1000LB01410  
 1000LB01411  
 1000LB01414

**Keypads and labels (continued)**



## Keypads and labels (continued)

**Note:**

Label **A** is a blank combined serial number and status label. Transfer information from old label. This label should be used in conjunction with the clear protective cover from the universal label set.

The picture above shows the label set that is available as a separate item from the standard Alaris® Syringe Pump label sets.

| Item | Description                                 | Part Number                                  |
|------|---------------------------------------------|----------------------------------------------|
| ABC  | ASENA Label Combined Serial Number & Status | 1000LB00102                                  |
|      |                                             | Use in conjunction with universal label set. |

# ***Appendix A***

## ***Specifications***

### **In this appendix**

|                                      |           |
|--------------------------------------|-----------|
| <b>Infusion</b>                      | <b>63</b> |
| <b>Electrical</b>                    | <b>64</b> |
| <b>Physical</b>                      | <b>64</b> |
| <b>Environmental</b>                 | <b>64</b> |
| <b>Electromagnetic Compatibility</b> | <b>65</b> |

## Specifications

To be used for reference only, for more detailed specifications refer to relevant Directions For Use. These specifications refer to all models covered by this manual unless otherwise stated.

### Infusion

| Infusion rate (ml/h) | Models <b>GH/CC/TIVA/PK</b> | Model <b>GS</b> |
|----------------------|-----------------------------|-----------------|
| <b>5ml syringes</b>  | 0.1 - 150                   | 0.1 - 150       |
| <b>10ml syringes</b> | 0.1 - 300                   | 0.1 - 200       |
| <b>20ml syringes</b> | 0.1 - 600                   | 0.1 - 200       |
| <b>30ml syringes</b> | 0.1 - 900                   | 0.1 - 200       |
| <b>50ml syringes</b> | 0.1 - 1200                  | 0.1 - 200       |

**Volume Infused** 0.0 - 9990ml

| Bolus rate (ml/h)    | Models <b>GH/CC/PK</b> | Model <b>GS</b> | Model <b>TIVA</b>          |
|----------------------|------------------------|-----------------|----------------------------|
| <b>5ml syringes</b>  | 10 - 150               | 10 - 150        | 150                        |
| <b>10ml syringes</b> | 10 - 300               | 10 - 300        | 150, 300                   |
| <b>20ml syringes</b> | 10 - 600               | 10 - 500        | 150, 300 or 600            |
| <b>30ml syringes</b> | 10 - 900               | 10 - 500        | 150, 300, 600 or 900       |
| <b>50ml syringes</b> | 10 - 1200              | 10 - 500        | 150, 300, 600, 900 or 1200 |

**Bolus volume limit** 0.5 (0.1ml - v2.3.x & above or v1.9.x); - 25.0ml

**Purge rate (ml/h)** 100 - 500

**Purge volume** 0.5 - 5ml



**During PURGE/BOLUS the pressure limit alarms are temporarily increased to the maximum settings level.**

**Keep Vein Open (KVO) rate**

0.1ml/h - 2.5ml/h

**End Of Syringe rate**

Stop, KVO (0.1ml/h to 2.5ml/h), or set rate if lower than KVO.

**Volume To Be Infused (VTBI)**

0.1ml - 1 00ml (0.1ml - 1000ml - v2.3.x & above or v1.9.x),, 1min - 24h

(Models CC/GH)

**VTBI done rate (Models CC/GH)**

Stop, KVO (0.1ml/h to 2.5ml/h), set rate if lower than KVO or continue at set rate.

**Near End Of Infusion alarm**

1min - 15min to end of infusion, or 10% of syringe volume, whichever is smaller.

**End Of Infusion (EOI) alarm**

0.1% - 5% of syringe volume.

**Pumping pressure limit**

L-1 to L-10 (100mmHg to 1000mmHg approx)

L-0 (50mmHg approx)

**Occlusion accuracy with pressure set (% of full scale) (Model CC):**

| Pressure (mmHg)                 | <b>0</b> | <b>25</b> | <b>500</b> | <b>1000</b> |
|---------------------------------|----------|-----------|------------|-------------|
| <b>At 23°C</b>                  | ±2%      | ±4%       | ±5%        | ±6%         |
| <b>Between 5°C and 40°C</b> ±4% | ±7%      | ±7%       | ±10%       |             |

**System accuracy**

Volumetric Mean ± 2% (nominal)

## Specifications

### Electrical

|                          |                                                                            |
|--------------------------|----------------------------------------------------------------------------|
| <b>Battery type</b>      | NiMH sealed rechargeable 7.2V/2.7Ah.                                       |
| <b>Battery life</b>      | Typically 4h from fully charged @ 5.0ml/h & 20°C under nominal conditions. |
| <b>Battery charging</b>  | 2.5h from fully discharged to 90% charged.                                 |
| <b>Fuse type</b>         | 2 x T-1.25A, slow blow.                                                    |
| <b>AC power supply</b>   | 115/230VAC, 50/60Hz, 20VA (nominal).                                       |
| <b>Electrical safety</b> | Class I Type CF.                                                           |

### Physical

|                   |                                                            |          |          |          |
|-------------------|------------------------------------------------------------|----------|----------|----------|
| <b>Weight</b>     | 2.7 kg (excluding power cable)                             |          |          |          |
| <b>Dimensions</b> | <b>Models</b>                                              | <b>W</b> | <b>H</b> | <b>D</b> |
|                   | <b>GS/GH/PK</b>                                            | 310mm    | 121mm    | 200mm    |
|                   | <b>CC</b>                                                  | 335mm    | 121mm    | 200mm    |
|                   | <b>TIVA</b>                                                | 310mm    | 121mm    | 200mm    |
| <b>Latex</b>      | The Alaris® Syringe Pump range does not contain any Latex. |          |          |          |

### Environmental

|                   |      |
|-------------------|------|
| <b>IPX Rating</b> | IPX1 |
|-------------------|------|

#### Operating limits

|                             |                   |
|-----------------------------|-------------------|
| <b>Temperature</b>          | +5°C to +40°C     |
| <b>Relative humidity</b>    | 20% to 90%        |
| <b>Atmospheric pressure</b> | 700hPa to 1060hPa |

#### Transport limits

|                             |                   |
|-----------------------------|-------------------|
| <b>Temperature</b>          | -30°C to +50°C    |
| <b>Relative humidity</b>    | 10% to 95%        |
| <b>Atmospheric pressure</b> | 500hPa to 1060hPa |

## Electromagnetic Compatibility

### Warning:

- The use of any accessory, transducer, or cable with the Alaris® Syringe Pump other than those specified may result in increased emissions or decreased immunity of the pump.
- The Alaris® Syringe Pump should not be used adjacent to or stacked with other equipment, however if adjacent or stacked use is necessary, the Alaris® Syringe Pump should be observed to verify normal operation in the configuration in which it will be used.

### Caution:

- The Alaris® Syringe Pump is a CISPR 11 Group 1 Class A Medical Equipment System and intended for use by healthcare professionals only.
- Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed, put into service and used according to the EMC information provided in the accompanying documents.
- Portable and Mobile RF communications can affect Medical Electrical Equipment.
- Operating the pump near equipment which radiates high energy radio frequencies (electro surgical or cauterizing equipment, portable radios, cellular telephones, etc.) may cause false alarm conditions. If this happens, reposition the pump away from the source of interference or turn off the pump and manually regulate the flow.

| Guidance and Manufacturer's Declaration – Electromagnetic Emissions                                                                                                                                                       |            |                                                                                                                                                                                                           |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>The Alaris® Syringe Pump is intended for use in the electromagnetic environment specified below.</p> <p>The customer or the user of the Alaris® Syringe Pump should assure that it is used in such an environment.</p> |            |                                                                                                                                                                                                           |
| Emissions Test                                                                                                                                                                                                            | Compliance | Electromagnetic Environment – Guidance                                                                                                                                                                    |
| <b>CISPR 11</b><br>RF Emissions                                                                                                                                                                                           | Group 1    | The pump uses RF energy only for its internal function in the normal product offering. Therefore, its RF emissions are very low and are not likely to cause any interface in nearby electronic equipment. |
| <b>CISPR 11</b><br>RF Emissions                                                                                                                                                                                           | Class A    | The pump is suitable for use in all establishments, other than domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.  |
| <b>EN 61000-3-2</b><br>Harmonic Emissions                                                                                                                                                                                 | Class A    |                                                                                                                                                                                                           |
| <b>EN 61000-3-3</b><br>Voltage Fluctuations,<br>Flicker Emissions                                                                                                                                                         | Complies   |                                                                                                                                                                                                           |

## Electromagnetic Compatibility (continued)

### Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The Alaris® Syringe Pump is intended for use in the electromagnetic environment specified below.

The customer or the user of Alaris® Syringe Pump should assure that it is used in such an environment.

| Immunity Test                                                                              | EN 60601-1-2 Test Level                                      | Compliance Level                                   | Electromagnetic Environment – Guidance                                                                                                                                                                                                                                                                                                          |
|--------------------------------------------------------------------------------------------|--------------------------------------------------------------|----------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>EN 61000-4-2</b> Electro-Static Discharge (ESD)                                         | ±6 kV contact<br>±8 kV air                                   | ±8 kV contact (Note 2)<br>±15 kV air (Note 2)      | Floors should be wood, concrete, or ceramic tile.<br>If floors are covered with synthetic material, the relative humidity should be at least 30 %.                                                                                                                                                                                              |
| <b>EN 61000-4-4</b><br>Electrical Fast Transient, Burst (EFT) (Note 3)                     | ±2 kV for power supply lines<br>±1 kV for input/output lines | ±2 kV for power supply lines<br>N/A (Note 4)       | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                                                                             |
| <b>EN 61000-4-5</b><br>Power Line Surge (Note 3)                                           | ±1 kV Line(s) to Line(s)<br>±2 kV Line(s) to Earth           | ±1 kV Line(s) to Line(s)<br>±2 kV Line(s) to Earth | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                                                                             |
| <b>EN 61000-4-8</b> Power Frequency Magnetic Field (50/60 Hz)                              | 3 A/m                                                        | 400 A/m 50 Hz (Note 2)                             | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                                                                                       |
| <b>EN 61000-4-11</b><br>Voltage Dips, Short Interruptions, and Voltage Variations (Note 3) | <5 % UT (Note 1)<br>(>95 % dip in UT)<br>for 0.5 cycle       | <5 % UT<br>(>95 % dip in UT)<br>for 0.5 cycle      | Mains power quality should be that of a typical commercial or hospital environment.<br><br>If the user of the pump requires continued operation during power mains interruptions, it is recommended that the pump be powered from an uninterruptible power supply or a battery.<br><br>The pump does employ an internal short duration battery. |
|                                                                                            | 40 % UT<br>(60 % dip in UT)<br>for 5 cycles                  | 40 % UT<br>(60 % dip in UT)<br>for 5 cycles        |                                                                                                                                                                                                                                                                                                                                                 |
|                                                                                            | 70 % UT<br>(30 % dip in UT)<br>for 25 cycles                 | 70 % UT<br>(30 % dip in UT)<br>for 25 cycles       |                                                                                                                                                                                                                                                                                                                                                 |
|                                                                                            | <5 % UT<br>(>95 % dip in UT)<br>for 5 sec                    | <5 % UT<br>(>95 % dip in UT)<br>for 5 sec          |                                                                                                                                                                                                                                                                                                                                                 |

Note 1— $U_T$  is the AC mains voltage prior to application of the test level.

Note 2—Compliance levels raised by EN 60601-2-24.

Note 3—Performed at the Minimum and Maximum Rated Input Voltage.

Note 4—Cardinal Health recommends using signal cables of less than 3 meters in length and this requirement is applicable only if signal cables are 3 meters or more in length. (EN 60601-1-2:2002, Clause 36.202.4)

## Electromagnetic Compatibility (continued)

### Guidance and Manufacturer's Declaration—Electromagnetic Immunity LIFE SUPPORT Equipment

The Alaris® Syringe Pump is intended for use in the electromagnetic environment specified below.  
The customer or the user of the Alaris® Syringe Pump should ensure that it is used in such an environment.

| Immunity Test                       | EN 60601-1-2<br>Test Level      | Compliance<br>Level  | Electromagnetic Environment – Guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|-------------------------------------|---------------------------------|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>EN 61000-4-6</b><br>Conducted RF | 3 V rms<br>150 kHz to 80<br>MHz | 10 V rms<br>(Note 3) | Portable and mobile RF communications equipment should be used no closer to any part of the pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><br><b>Recommended Separation Distance</b><br><br>$d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>EN 61000-4-3</b><br>Radiated RF  | 3 V/m<br>80 MHz to 2.5 GHz      | 10 V/m<br>(Note 3)   | $d = \left[ \frac{12}{V_2} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{12}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 2.5 \text{ GHz}$ $d = \left[ \frac{23}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).<sup>a</sup></p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, b should be less than the compliance level in each frequency range. c</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |

Note 1—At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Note 3—Compliance levels raised by EN 60601-2-24.

<sup>a</sup> The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.

<sup>b</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the pump is used exceeds the applicable RF compliance level above, the pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the pump.

<sup>c</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

## Electromagnetic Compatibility (continued)

### Recommended Separation Distances for LIFE SUPPORT Equipment between portable and mobile RF communications equipment and the Alaris® Syringe Pump

The Alaris® Syringe Pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled.

The user of the Alaris® Syringe Pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Alaris® Syringe Pump as recommended below, according to the maximum output power of the communications equipment.

| Rated Maximum Output Power of Transmitter<br>W | Separation Distance According to Frequency of Transmitter<br>m                              |                                                                                      |                                                                      |                                                                       |
|------------------------------------------------|---------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|----------------------------------------------------------------------|-----------------------------------------------------------------------|
|                                                | 150 kHz to 80 MHz<br>Outside ISM bands 3.5<br>$d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ | 150 kHz to 80 MHz<br>In ISM bands 12<br>$d = \left[ \frac{12}{V_2} \right] \sqrt{P}$ | 80 MHz to 800 MHz 12<br>$d = \left[ \frac{12}{E_1} \right] \sqrt{P}$ | 800 MHz to 2.5 GHz 23<br>$d = \left[ \frac{23}{E_1} \right] \sqrt{P}$ |
| 0.01                                           | 0.03                                                                                        | 0.12                                                                                 | 0.12                                                                 | 0.23                                                                  |
| 0.1                                            | 0.11                                                                                        | 0.38                                                                                 | 0.38                                                                 | 0.73                                                                  |
| 1                                              | 0.35                                                                                        | 1.20                                                                                 | 1.20                                                                 | 2.30                                                                  |
| 10                                             | 1.11                                                                                        | 3.80                                                                                 | 3.80                                                                 | 7.28                                                                  |
| 100                                            | 3.50                                                                                        | 12.00                                                                                | 12.00                                                                | 23.00                                                                 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1—At 80 MHz and 800 MHz, the separation distance for the higher frequency range apply.

Note 2—The ISM (Industrial, Scientific, and Medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

Note 3—An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

Note 4—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

# Appendix B

## Disposal

### In this appendix

|                 |    |
|-----------------|----|
| Disposal        | 70 |
| Battery Removal | 70 |

## Disposal

### Information on Disposal for Users of Waste Electrical & Electronic Equipment

This  symbol on the product and/or accompanying documents means that used electrical and electronic products should not be mixed with municipal waste.

If you wish to discard electrical and electronic equipment, please contact your Cardinal Health affiliate office or distributor for further information.

Disposing of this product correctly will help to save valuable resources and prevent any potential negative effects on human health and the environment which could otherwise arise from inappropriate waste handling.


### Information on Disposal in Countries outside the European Union

This symbol is only valid in the European Union. The product should be disposed of taking environmental factors into consideration. To ensure no risk or hazard, remove the internal rechargeable battery and the Nickel Metal Hydride battery from the control board and dispose of as outlined by the local country regulations. All other components can be safely disposed of as per local regulations.

## Battery Removal

### Remove the Main Battery

Remove the two case screws in battery cover, remove cover and battery.




### Remove the Battery on Control PCB

Remove the Control PCB from the pump, see 'Spare Parts Replacement Procedures'.

Remove the protective plastic cover from the Control PCB.  
Desolder battery from the Control PCB.

#### Control PCB Reverse side



## *Spare Parts Listing*

### **In this appendix**

|                                   |           |
|-----------------------------------|-----------|
| <b>Electrical parts listing</b>   | <b>72</b> |
| <b>Front case parts listing</b>   | <b>72</b> |
| <b>Rear case parts listing</b>    | <b>73</b> |
| <b>Keypads and labels</b>         | <b>73</b> |
| <b>Transmission parts listing</b> | <b>74</b> |
| <b>Software</b>                   | <b>74</b> |
| <b>Test equipment</b>             | <b>74</b> |

## Electrical Parts Listing

| Part Number | DESCRIPTION                                |
|-------------|--------------------------------------------|
| 1000EL00222 | ASENA SP, FUSE, T-1.25A SLOW BLOW, MAINS   |
| 1000SP00189 | ASENA SP, KIT, CHASSIS PCB                 |
| 1000SP00497 | ASENA TIVA, CONTROL PCB (SOFTWARE V2.2.1)  |
| 1000SP01118 | ASENA GS, DISPLAY PCB                      |
| 1000SP01119 | ASENA GH/CC/TIVA, DISPLAY PCB              |
| 1000SP00495 | ASENA GH, CONTROL PCB (SOFTWARE V2.2.1)    |
| 1000SP01122 | ASENA SP, ASSY, BATTERY                    |
| 1000SP01124 | ASENA SP, KIT, MAINS INLET                 |
| 1000SP01125 | ASENA SP, ASSY, SYRINGE SIZE POTENTIOMETER |
| 1000SP01130 | ASENA SP, KIT, SPEAKER                     |
| 1000SP00494 | ASENA GS, CONTROL PCB (SOFTWARE V2.2.1)    |
| 1000SP01155 | ASENA CC, KIT, PRESSURE TRANSDUCER         |
| 1000SP01160 | ASENA SP, KIT, RS232 (PCB & FIXINGS)       |
| 1000SP00496 | ASENA CC, CONTROL PCB (SOFTWARE V2.2.1)    |
| 8000EL00019 | ASENA SP, PLUNGER DETECT PCB               |
| 8000EL00022 | ASENA SP, CARRIAGE PCB                     |
| 8000EL00063 | ASENA SP, ASSY, POWER SUPPLY UNIT (PSU)    |

## Front Case Parts Listing

| Part Number | DESCRIPTION                                 |
|-------------|---------------------------------------------|
| 1000SP01153 | ASENA CC, KIT, FRONT CASE                   |
| 1000SP00478 | ASENA GS, KIT, FRONT CASE                   |
| 1000SP00479 | ASENA GH, KIT, FRONT CASE                   |
| 1000SP00480 | ASENA TIVA , KIT, FRONT CASE                |
| 1000SP01204 | ASENA PK , KIT, FRONT CASE                  |
| 1000SP01167 | ASENA SP, KIT, SPARE FOOT                   |
| 1000SP01137 | ASENA GH/CC/TIVA, KIT, DISPLAY PROTECTION   |
| 1000SP01123 | ASENA SP, KIT, SYRINGE & FLANGE CLAMPS      |
| 1000SP01116 | ASENA SP, KIT, SYRINGE SIZER                |
| 1000SP00122 | ASENA GS, ASSY, WINDOW ADAPTOR DISPLAY      |
| 1000SP00187 | ASENA GS, ASSY, DISPLAY INSULATOR           |
| 1000SP00188 | ASENA GH/CC/TIVA, ASSY, DISPLAY INSULATOR   |
| 0000ME00423 | PAD SELF ADHESIVE DOUBLE SIDED 12X12MM      |
| 1000ME00261 | ASENA SP, ASSY, BRACKET DISPLAY MOUNTING    |
| 1000ME01301 | ASENA SP, ASSY, GASKET DISPLAY              |
| 1000ME00450 | ASENA CC, ASSY, TOP DISC HOLDER             |
| 1000ME01318 | ASENA CC, ASSY, FLAG DISC HOLDER            |
| 1000ME01319 | ASENA CC, ASSY, GASKET MYLAR                |
| 1000ME00311 | ASENA SP, CASE SEALING CORD (1M)            |
| 1000SP00466 | ASENA SP, KIT, FIXINGS (SCREWS,WASHERS,ETC) |

## Rear Case Parts Listing

| Part Number | DESCRIPTION                                          |
|-------------|------------------------------------------------------|
| 1000SP01115 | ASENA GS/GH/TIVA, KIT, REAR CASE                     |
| 1000SP01154 | ASENA CC, KIT, REAR CASE                             |
| 1000SP01121 | ASENA SP, BATTERY COVER/HANDLE                       |
| 1000SP00593 | ASENA SP, KIT, SPARE ADHESIVE FOOT RIVET REPLACEMENT |
| 1000SP00595 | ASENA SP, KIT, SPARE ADHESIVE FOOT                   |
| 1000SP01114 | ASENA SP, KIT, RAIL CAM                              |
| 1000SP00115 | ASENA SP, ASSY, POLE CLAMP                           |
| 1000ME01213 | ASENA SP, TUBE RESTRAINT BLANK                       |
| 1000ME01214 | ASENA SP, TUBE RESTRAINT RS232                       |
| 1000ME01303 | MAGNET IR DETECT                                     |
| 1000ME01306 | ASENA SP, ASSY, PSU INSULATOR                        |
| 1000ME01317 | ASENA CC, ASSY, PLUG BLANKING TRANSDUCER             |
| 1000SP00466 | ASENA SP, KIT, FIXINGS (SCREWS,WASHERS,ETC)          |
| 1000SP00467 | ASENA SP/GW, KIT, PE STUD                            |
| 1000SP00468 | ASENA SP/GW, KIT, RS232 CONNECTOR                    |
| 1000SP00589 | SPARE KIT POLE CLAMP ARM                             |
| 0000ME00379 | SPRING POLE CLAMP                                    |
| 0000ME00369 | BEARING BALL 5MM SS AISI316 GRADE 100                |
| 1000ME01244 | CLAMP POLE PIVOT (SHOULDER SCREW)                    |

## Keypads and Labels

| Part Number | DESCRIPTION                              |
|-------------|------------------------------------------|
| 1000SP01126 | ASENA GH/CC, KEY, SWITCH OVERLAY KEYPAD  |
| 1000SP01127 | ASENA GS, KEY, SWITCH OVERLAY KEYPAD     |
| 1000SP00403 | ASENA TIVA, KEY, SWITCH OVERLAY KEYPAD   |
| 1000SP01203 | ASENA PK, KEY, SWITCH OVERLAY KEYPAD     |
| 1000LB00431 | ASENA SP, LBL, LABEL CHASSIS SCREW COVER |
| 1000LB01408 | LABEL SET CC SYRINGE PUMP                |
| 1000LB01409 | LABEL SET GH SYRINGE PUMP                |
| 1000LB01410 | LABEL SET GS SYRINGE PUMP                |
| 1000LB01411 | LABEL SET TIVA SYRINGE PUMP              |
| 1000LB01414 | LABEL SET PK SYRINGE PUMP                |

## Transmission Parts Listings

| Part Number | DESCRIPTION                            |
|-------------|----------------------------------------|
| 1000ME01218 | ASENA SP, GRIPPER BOTTOM               |
| 1000ME01219 | ASENA SP, GRIPPER TOP                  |
| 1000SP00230 | ASENA SP, KIT, CARRIAGE BUFFER         |
| 1000SP01109 | ASENA SP, KIT, STEPPER MOTOR           |
| 1000SP01112 | ASENA SP, KIT, CHASSIS ASSEMBLY        |
| 1000SP01113 | ASENA SP, KIT, PLUNGER ASSEMBLY        |
| 1000SP01136 | ASENA SP, KIT, CHASSIS ENHANCEMENT     |
| 1000SP01110 | ASENA SP, KIT, MOTOR PLATE             |
| 1000ME01451 | ASENA SP, ASSY, LINEAR POTENTIOMETER   |
| 1000SP01107 | SPARE CARRIAGE ASENA                   |
| 1000ME01198 | GEAR DECLUTCH                          |
| 1000SP00408 | SP KIT MOTOR PLATE STRAIN BEAM SUPPORT |

## Software

| Part Number | DESCRIPTION                                            |
|-------------|--------------------------------------------------------|
| 1000SP01223 | ASENA CC, SOFT, SOFTWARE CD V1.5.10 (MK1 & 2)          |
| 1000SP01221 | ASENA GS, SOFT, SOFTWARE CD V1.5.10 (MK1 & 2)          |
| 1000SP01222 | ASENA GH, SOFT, SOFTWARE CD V1.5.10 (MK1 & 2)          |
| 1000SP01224 | ASENA TIVA, SOFT, SOFTWARE CD V1.6.2 (MK1 & 2)         |
| 1000SP01270 | ASENA SYRINGE PUMP, SOFT, SOFTWARE CD V1.9.4 (MK1 & 2) |
| 1000SP01227 | ASENA CC, SOFT, SOFTWARE CD V2.0.0 (MK3)               |
| 1000SP01225 | ASENA GS, SOFT, SOFTWARE CD V2.0.0 (MK3)               |
| 1000SP01226 | ASENA GH, SOFT, SOFTWARE CD V2.0.0 (MK3)               |
| 1000SP01228 | ASENA TIVA, SOFT, SOFTWARE CD V2.1.0 (MK3)             |
| 1000SP01267 | ASENA CC, SOFT, SOFTWARE CD V2.3.6 (MK3)               |
| 1000SP01276 | ASENA GS, SOFT, SOFTWARE CD V2.3.6 (MK3)               |
| 1000SP01268 | ASENA GH, SOFT, SOFTWARE CD V2.3.6 (MK3)               |
| 1000SP01269 | ASENA TIVA, SOFT, SOFTWARE CD V2.3.6 (MK3)             |
| 1000SP01202 | ASENA PK, SOFT, SOFTWARE CD V3.2.12 (MK3)              |

## Test Equipment

| Part Number | DESCRIPTION                                 |
|-------------|---------------------------------------------|
| 0000JG00047 | ASENA SP, TEST, CASE CRADLE JIG             |
| 0000JG00053 | ASENA SP, TEST, PE STUD SOCKET              |
| 0000TG00080 | ASENA SP & P SERIES, TEST, RATE GAUGE       |
| 1000TG00095 | LINEAR SIZING SPACER BOM                    |
| 1000SP00065 | ASENA SP & P SERIES, KIT, CALIBRATION TOOLS |
| 1000SP00172 | ASENA SP, KIT, IRDA PORT CABLE & HEADER PCB |
| 1000SP00209 | ASENA SP, KIT, EVENT LOG DOWNLOAD UTILITY   |
| 0000JG00014 | ASENA SP & P SERIES, TEST, PLUNGER PROTECT  |
| 1000SP00336 | ASENA, ASSY, RS232 CABLE                    |
| 1000SP00481 | ASENA SP, TEST, SYRINGE CLAMP JIG           |
| 1000ME01466 | POLE CLAMP SNAKE EYE DRIVER                 |
| 0000TG00200 | DIGITAL OCCLUSION TEST GEAR                 |

# ***Appendix D***

## ***Fitting & Replacement Guidelines***

### **In this appendix**

|                                     |           |
|-------------------------------------|-----------|
| <b>General assembly information</b> | <b>76</b> |
| <b>Torque guide</b>                 | <b>76</b> |

## General assembly information

1. A wide range of self-tapping fasteners are available.
2. PT screws are for plastic, self-tapping applications.
3. TAPTITE screws are for metal self thread-forming applications. These can be recognised by a triangular cross-section on the end.
4. Almost all fasteners on the Alaris® Syringe Pumps are self-tapping and have the potential to be over-tightened (over-torqued).
5. The force required to create a thread for the first time is more than when reassembling a previously made joint.
6. Always use the correct torque level when first making an assembly stage.
7. Take care with the torque applied when re-assembling parts. Less torque is required, so a hand tool may be more appropriate.
8. In many situations a stripped thread will require replacement of the failed component.
9. The head patterns of the fasteners are of the following types:
  - ◆ Pozi Number 1 (smaller X head)
  - ◆ Pozi Number 2 (larger X head)
  - ◆ Torx Number T8 (Small star profile, used typically on countersunk parts with smaller heads - Backplate / Mains Inlet / Carriage PCB screw).
  - ◆ Torx Number T10 (Medium star profile, used on the majority of Alaris® Syringe Pump Torx fasteners)
  - ◆ Torx Number T20 (Larger star shape, typically for case securing screws)
  - ◆ M3 (Hex head with 5.5mm across flats (AF) drivers)
  - ◆ M4 nuts (Hex head with 7mm across flats (AF) drivers)
10. Always select the correct tool and bit pattern for the fastener.

## Torque guide

1. When selecting a torque for a servicing activity, be aware that refastening will require less torque than the initial manufacture.
2. Use this information as a guide to the 'do not exceed' torque levels when servicing the equipment. When servicing it is recommended that torque is applied gradually until the component is secure. In any process do not exceed the stated levels.
3. If a torque driver is available for servicing this will help control the applied torque. Otherwise, be aware that excess force may cause the component to fail.

### Plunger Drive Assembly:

| Stage Description               | Component Description                  | Qty | Established Process Torque |
|---------------------------------|----------------------------------------|-----|----------------------------|
| Intermediate Tube Bearing Plate | Screw - PT K30x8 Pan Hd Torx (T10)     | 2   | 50 cNm                     |
| Fixing Gripper Gears            | Screw - PT K22x12 Pan Hd (T6)          | 2   | 50 cNm                     |
| Screw on the Backplate Assembly | Screw - PT K30x12 Csk Torx (T8) Rogard | 3   | 40 cNm                     |

### Main Chassis Assembly:

| Stage Description                      | Component Description              | Qty | Established Process Torque |
|----------------------------------------|------------------------------------|-----|----------------------------|
| Mount motor onto motor plate           | Screw - M3x12 Pan Hd Torx (T10)    | 3   | 60 cNm                     |
| Mount motor plate to chassis assembly  | Screw - Taptite M4x10 Csk Pozi     | 3   | 1.0 Nm                     |
| Attach drive belt and leadscrew pulley | Nut - M4 Full                      | 1   | 40 cNm                     |
| Secure Carriage PCB to carriage        | Screw - PT K30x6 Csk Pozi          | 1   | 30 cNm                     |
| Secure carriage plate to carriage      | Screw - PT K30x6 Pan Hd Torx (T10) | 1   | 40 cNm                     |
| Secure bearing block to chassis        | Screw - Taptite M4x10 Csk Pozi     | 3   | 1.0 Nm                     |

## Torque guide (continued)

### Front Case Assembly:

| Stage Description                                          | Component Description               | Qty | Established Process Torque |
|------------------------------------------------------------|-------------------------------------|-----|----------------------------|
| Secure syringe sizing retainer                             | Screw - PT K30x8 Csk Pozi           | 2   | 40 cNm                     |
| Attach syringe clamp                                       | Screw M3 x 8 Pan Hd Torx (T10)      | 1   | 50 cNm                     |
| Secure display / mounting brackets to front case           | Screw PT K30x12 Csk Pozi            | 4   | 50 cNm                     |
| Secure syringe flange clamp to chassis / bearing block     | Screw - PT K30x14 Pan Hd Torx (T10) | 2   | 70 cNm                     |
| Secure chassis to front case                               | Screw M3x8 Taptite Csk Pozi         | 1   | 50 cNm                     |
| Secure plunger drive to carriage                           | Screw - PT K30x6 Pan Hd Torx (T10)  | 2   | 30 cNm                     |
| Secure RS232 option                                        | Spacer 8mm Hex Br/Ni PI M3x6mm      | 4   | hand tight                 |
|                                                            | Nut M3 St. St A2                    | 4   | 40 cNm                     |
|                                                            | Screw M3x6 Pan Hd Pozi Z+C          | 4   | 40 cNm                     |
|                                                            | Screw - PT K30x6 Pan Hd Torx (T10)  | 3   | 30 cNm                     |
| Secure Control board                                       | Screw - M3x8 Taptite Pan Hd         | 2   | 50 cNm                     |
| Secure Chassis PCB to chassis                              | Screw - PT K30x12 Pan Hd Torx (T10) | 2   | 70 cNm                     |
| Secure Model CC pressure transducer assembly to front case |                                     |     |                            |

### Assembly Rear Case:

| Stage Description                                         | Component Description                  | Qty | Established Process Torque |
|-----------------------------------------------------------|----------------------------------------|-----|----------------------------|
| Attach pole clamp to rear case                            | Screw - M3x8 Pan Hd Torx (T10)         | 3   | 70 cNm                     |
| Attach pole clamp arm to pole v clamp                     | Pivot screw                            | 1   | 2 Nm                       |
| Secure rail clamp to cam                                  | Screw - PT K30x10 Csk Torx (T8) Rogard | 1   | 70 cNm                     |
| Attach camera lever                                       | Screw PT K30x8 Pan Hd Torx (T10)       | 1   | 60 cNm                     |
| Secure mains inlet assy to retainer                       | Screw - PT K30x12 Csk Torx (T8) Rogard | 2   | 40 cNm                     |
| Secure PE stud                                            | M6 Nut                                 | 2   | hand tight                 |
| Secure PSU to case                                        | Screw - PT K30x8 Pan Hd Torx (T10)     | 3   | 40 cNm                     |
| Secure earth lead to PSU metal frame                      | Screw M3x6 Pan Hd Pozi                 | 1   | hand tight                 |
| Fit tube restraint / RS232 cover                          | Screw PT K30x8 Csk                     | 2   | 40 cNm                     |
| Fit RS232 male / female connectors                        | Jack Socket RS232 P8000                | 1   | hand tight                 |
| Secure Model CC pressure transducer assembly to rear case | Screw - 4gx1/2" S/T B Zn Clr Pz1       | 1   | 60 cNm                     |

### Final Assembly:

| Stage Description              | Component Description             | Qty | Established Process Torque |
|--------------------------------|-----------------------------------|-----|----------------------------|
| Secure front case to rear case | Screw PT K40x12 Pan Hd Torx (T20) | 6   | 70 cNm                     |
| Secure battery cover / handle  | Screw PT K40x12 Pan Hd Torx (T20) | 2   | 70 cNm                     |

### Pressure Transducer Assembly for Model CC only:

| Stage Description                 | Component Description            | Qty | Established Process Torque |
|-----------------------------------|----------------------------------|-----|----------------------------|
| Secure disk holder centre to top  | Screw PT K30x8 Pan Hd Torx (T10) | 1   | 50 cNm                     |
| Secure disk holder base to centre | Screw PT K30x8 Pan Hd Torx (T10) | 4   | 50 cNm                     |

# ***Appendix E***

## ***Configuration & Drug Protocol Records***

### **In this appendix**

|                                                            |           |
|------------------------------------------------------------|-----------|
| <b>Alaris® GS Syringe Pump Configured Options Record</b>   | <b>79</b> |
| <b>Alaris® GH Syringe Pump Configured Options Record</b>   | <b>80</b> |
| <b>Alaris® CC Syringe Pump Configured Options Record</b>   | <b>81</b> |
| <b>Alaris® TIVA Syringe Pump Configured Options Record</b> | <b>82</b> |
| <b>Alaris® CC Syringe Pump Drug Protocol Setup</b>         | <b>83</b> |
| <b>Alaris® TIVA Syringe Pump Drug Protocol Setup</b>       | <b>84</b> |

## Alaris® GS Syringe Pump Configured Options Record

### General Options

Enter the pump-specific information for your records on a copy of this page.

| Option             | Default           |                          | Range                       | Setting |
|--------------------|-------------------|--------------------------|-----------------------------|---------|
| Software Version   | 1.5.10 & 2.0.0    | 1.9.x<br>2.3.x and above |                             |         |
| NURSE CALL FITTED  | Disabled          | Disabled                 | Enabled/Disabled            |         |
| NURSE CALL INVERT  | Disabled          | Disabled                 | Enabled/Disabled            |         |
| RS232 SELECTED     | Disabled          | Disabled                 | Enabled/Disabled            |         |
| NEOI WARNING       | 1min              | 5mins                    | 1min - 15mins               |         |
| EOI POINT          | 1.0%              | 1.0%                     | 0.1% - 5% of syringe volume |         |
| KVO AT EOI         | Enabled           | Enabled                  | Enabled/Disabled            |         |
| KVO RATE           | 1.0ml/h           | 1.0ml/h                  | 0.1ml/h - 2.5ml/h           |         |
| BACK OFF           |                   | Enabled                  | Enabled/Disabled            |         |
| AUTO SAVE          | Enabled           | Enabled                  | Enabled/Disabled            |         |
| RATE LOCK          | Disabled          | Disabled                 | Enabled/Disabled            |         |
| QUIET MODE         | Disabled          | Disabled                 | Enabled/Disabled            |         |
| AC FAIL            | Enabled           | Enabled                  | Enabled/Disabled            |         |
| PRESSURE DISPLAY   | Disabled          | Enabled                  | Enabled/Disabled            |         |
| PRESSURE DEFAULT   | L3                | L3                       | L0 - L10(50mmHg -1000mmHg)  |         |
| CAP RATE           | Max infusion rate | 200ml/h                  | 1.0ml/h - 200ml/h           |         |
| PURGE RATE         | 200ml/h           | 200ml/h                  | 100ml/h - 500ml/h           |         |
| PURGE VOLUME LIMIT | 2.0ml             | 2.0ml                    | 0.5ml - 5.0ml               |         |
| PURGE SYRINGE      |                   | Disabled                 | Enabled/Disabled            |         |
| BOLUS              | Enabled           | Enabled                  | Enabled/Disabled            |         |
| DEFAULT BOLUS      | Max bolus rate    | 500ml/h                  | 10ml/h - 500ml/h            |         |
| CAP BOLUS RATE     | Max bolus rate    | 500ml/h                  | 10ml/h - 500ml/h            |         |
| BOLUS VOL LIMIT    | 5.0ml             | 5.0ml                    | 0.5ml (0.1ml)* - 25.0ml     |         |
| MANUAL BOLUS       |                   | Disabled                 | Enabled/Disabled            |         |
| CALL BACK TIME     |                   | 2.0mins                  | 0.1mins - 15mins            |         |
| EVENT LOG DISPLAY  | Disabled          | Enabled                  | Enabled/Disabled            |         |
| BATTERY ICON       |                   | Enabled                  | Enabled/Disabled            |         |
| AUDIO VOLUME       | Medium            | Medium                   | Low, Medium, High           |         |
| AUTO NIGHT MODE    | Enabled           | Enabled                  | Enabled/Disabled            |         |

### Syringes Enabled

\* For software versions 1.9.x & 2.3.x and above

| Make | Size(s) | Make | Size(s) |
|------|---------|------|---------|
|      |         |      |         |
|      |         |      |         |
|      |         |      |         |
|      |         |      |         |

**Hospital Name**

**Serial No.**

**Software Version**

**Approved by**

**Configured by**

**Date**

**Date**

## Configuration & Drug Protocol Records

### Alaris® GH Syringe Pump Configured Options Record

#### General Options

Enter the pump-specific information for your records on a copy of this page.

| Option             | Default           |                          | Range                       | Setting |
|--------------------|-------------------|--------------------------|-----------------------------|---------|
| Software Version   | 1.5.10 & 2.0.0    | 1.9.x<br>2.3.x and above |                             |         |
| NURSE CALL FITTED  | Disabled          | Disabled                 | Enabled/Disabled            |         |
| NURSE CALL INVERT  | Disabled          | Disabled                 | Enabled/Disabled            |         |
| RS232 SELECTED     | Disabled          | Disabled                 | Enabled/Disabled            |         |
| NEOI WARNING       | 1min              | 5mins                    | 1min - 15mins               |         |
| EOI POINT          | 1.0%              | 1.0%                     | 0.1% - 5% of syringe volume |         |
| KVO AT EOI         | Enabled           | Enabled                  | Enabled/Disabled            |         |
| KVO RATE           | 1.0ml/h           | 1.0ml/h                  | 0.1ml/h - 2.5ml/h           |         |
| BACK OFF           | Disabled          | Enabled                  | Enabled/Disabled            |         |
| AUTO SAVE          | Enabled           | Enabled                  | Enabled/Disabled            |         |
| RATE LOCK          | Disabled          | Disabled                 | Enabled/Disabled            |         |
| QUIET MODE         | Disabled          | Disabled                 | Enabled/Disabled            |         |
| AC FAIL            | Enabled           | Enabled                  | Enabled/Disabled            |         |
| RATE TITRATION     | Disabled          | Disabled                 | Enabled/Disabled            |         |
| PRESSURE DISPLAY   | Disabled          | Enabled                  | Enabled/Disabled            |         |
| CAP PRESSURE       |                   | L10                      | L0 - L10(50mmHg -1000mmHg)  |         |
| PRESSURE DEFAULT   | L3                | L3                       | L0 - L10(50mmHg -1000mmHg)  |         |
| CAP RATE           | Max infusion rate | 1200ml/h                 | 1.0ml/h - 1200ml/h          |         |
| PURGE RATE         | 200ml/h           | 200ml/h                  | 100ml/h - 500ml/h           |         |
| PURGE VOLUME LIMIT | 2.0ml             | 2.0ml                    | 0.5ml - 5.0ml               |         |
| PURGE SYRINGE      |                   | Disabled                 | Enabled/Disabled            |         |
| BOLUS              | Enabled           | Enabled                  | Enabled/Disabled            |         |
| DEFAULT BOLUS      | Max bolus rate    | 500ml/h                  | 10ml/h - 1200ml/h           |         |
| CAP BOLUS RATE     | Max bolus rate    | 1200ml/h                 | 10ml/h - 1200ml/h           |         |
| BOLUS VOL LIMIT    | 5.0ml             | 5.0ml                    | 0.5ml (0.1ml)* - 25.0ml     |         |
| MANUAL BOLUS       |                   | Disabled                 | Enabled/Disabled            |         |
| CALL BACK TIME     |                   | 2mins                    | 0.1mins - 15mins            |         |
| VTBI CLEAR RATE    |                   | Disabled                 | Enabled/Disabled            |         |
| EVENT LOG DISPLAY  | Disabled          | Enabled                  | Enabled/Disabled            |         |
| BATTERY ICON       |                   | Enabled                  | Enabled/Disabled            |         |
| AUDIO VOLUME       | Medium            | Medium                   | Low, Medium, High           |         |
| AUTO NIGHT MODE    | Enabled           | Enabled                  | Enabled/Disabled            |         |

#### Syringes Enabled

| Make | Size(s) |
|------|---------|
|      |         |
|      |         |
|      |         |
|      |         |
|      |         |

#### Drug Names

\* For software versions 1.9.x & 2.3.x and above

|   |  |    |  |
|---|--|----|--|
| 1 |  | 7  |  |
| 2 |  | 8  |  |
| 3 |  | 9  |  |
| 4 |  | 10 |  |
| 5 |  | 11 |  |
| 6 |  | 12 |  |

Hospital Name

Serial No.

Software Version

Approved by

Configured by

Date

Date

## Configuration & Drug Protocol Records

### Alaris® CC Syringe Pump Configured Options Record

#### General Options

Enter the pump-specific information for your records on a copy of this page.

| Option             | Default           |                          | Range                       | Setting |
|--------------------|-------------------|--------------------------|-----------------------------|---------|
| Software Version   | 1.5.10 & 2.0.0    | 1.9.x<br>2.3.x and above |                             |         |
| NURSE CALL FITTED  | Disabled          | Disabled                 | Enabled/Disabled            |         |
| NURSE CALL INVERT  | Disabled          | Disabled                 | Enabled/Disabled            |         |
| RS232 SELECTED     | Disabled          | Disabled                 | Enabled/Disabled            |         |
| NEOI WARNING       | 1min              | 5mins                    | 1min - 15mins               |         |
| EOI POINT          | 1.0%              | 1.0%                     | 0.1% - 5% of syringe volume |         |
| KVO AT EOI         | Enabled           | Enabled                  | Enabled/Disabled            |         |
| KVO RATE           | 1.0ml/h           | 1.0ml/h                  | 0.1ml/h - 2.5ml/h           |         |
| BACK OFF           | Disabled          | Enabled                  | Enabled/Disabled            |         |
| AUTO SAVE          | Enabled           | Enabled                  | Enabled/Disabled            |         |
| RATE LOCK          | Disabled          | Disabled                 | Enabled/Disabled            |         |
| QUIET MODE         | Disabled          | Disabled                 | Enabled/Disabled            |         |
| AC FAIL            | Enabled           | Enabled                  | Enabled/Disabled            |         |
| RATE TITRATION     | Disabled          | Disabled                 | Enabled/Disabled            |         |
| PRESSURE DISPLAY   | Disabled          | Enabled                  | Enabled/Disabled            |         |
| AUTO PRESSURE      | Disabled          | Enabled                  | Enabled/Disabled            |         |
| AUTO SET PRESSURE  |                   | Disabled                 | Enabled/Disabled            |         |
| AUTO OFFSET        |                   | 30mmHg                   | 15mmHg - 100mmHg            |         |
| PRESSURE DEFAULT   | 300mmHg           | 300mmHg                  | 1mmHg - 1000mmHg            |         |
| MAX PRESSURE       | 1000mmHg          | 1000mmHg                 | 1mmHg - 1000mmHg            |         |
| WEIGHT             | 70Kg              | 1.00Kg                   | 0.01Kg - 250Kg              |         |
| CAP RATE           | Max infusion rate | 1200ml/h                 | 1.0ml/h - 1200ml/h          |         |
| PURGE RATE         | 200ml/h           | 200ml/h                  | 100ml/h - 500ml/h           |         |
| PURGE VOLUME LIMIT | 2.0ml             | 2.0ml                    | 0.5ml - 5.0ml               |         |
| PURGE SYRINGE      |                   | Disabled                 | Enabled/Disabled            |         |
| BOLUS              | Enabled           | Enabled                  | Enabled/Disabled            |         |
| DEFAULT BOLUS      | Max bolus rate    | 500ml/h                  | 10ml/h - 1200ml/h           |         |
| CAP BOLUS RATE     | Max bolus rate    | 1200ml/h                 | 10ml/h - 1200ml/h           |         |
| BOLUS VOL LIMIT    | 5.0ml             | 5.0ml                    | 0.5ml (0.1ml)* - 25.0ml     |         |
| MANUAL BOLUS       |                   | Disabled                 | Enabled/Disabled            |         |
| CALL BACK TIME     |                   | 2.0mins                  | 0.1mins - 15mins            |         |
| VTBI CLEAR RATE    |                   | Disabled                 | Enabled/Disabled            |         |
| EVENT LOG DISPLAY  | Disabled          | Enabled                  | Enabled/Disabled            |         |
| BATTERY ICON       |                   | Enabled                  | Enabled/Disabled            |         |
| AUDIO VOLUME       | Medium            | Medium                   | Low, medium, high           |         |
| AUTO NIGHT MODE    | Enabled           | Enabled                  | Enabled/Disabled            |         |

#### Units Enabled

\* For software versions 1.9.x & 2.3.x and above

|                                    |                                    |                                    |                                    |                                    |                                   |                                 |                                   |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|-----------------------------------|---------------------------------|-----------------------------------|
| <input type="checkbox"/> ng/min    | <input type="checkbox"/> µg/kg/min | <input type="checkbox"/> µg/24h    | <input type="checkbox"/> mg/kg/min | <input type="checkbox"/> mg/24h    | <input type="checkbox"/> g/24h    | <input type="checkbox"/> U/h    | <input type="checkbox"/> U/kg/24h |
| <input type="checkbox"/> ng/kg/min | <input type="checkbox"/> µg/h      | <input type="checkbox"/> µg/kg/24h | <input type="checkbox"/> mg/h      | <input type="checkbox"/> mg/kg/24h | <input type="checkbox"/> U/min    | <input type="checkbox"/> U/kg/h | <input type="checkbox"/> kU/24h   |
| <input type="checkbox"/> µg/min    | <input type="checkbox"/> µg/kg/h   | <input type="checkbox"/> mg/min    | <input type="checkbox"/> mg/kg/h   | <input type="checkbox"/> g/h       | <input type="checkbox"/> U/kg/min | <input type="checkbox"/> U/24h  | <input type="checkbox"/> mmol/h   |

#### Syringes Enabled

| Make | Size(s) | Make | Size(s) |
|------|---------|------|---------|
|      |         |      |         |
|      |         |      |         |
|      |         |      |         |

Hospital Name

Serial No.

Software Version

Approved by

Configured by

Date

Date

## Alaris® TIVA Syringe Pump Configured Options Record

### General Options

Enter the pump-specific information for your records on a copy of this page.

| Option             | Default       |                          | Range                       | Setting |
|--------------------|---------------|--------------------------|-----------------------------|---------|
| Software Version   | 1.6.2 & 2.1.0 | 1.9.x<br>2.3.x and above |                             |         |
| NURSE CALL FITTED  | Disabled      | Disabled                 | Enabled/Disabled            |         |
| NURSE CALL INVERT  | Disabled      | Disabled                 | Enabled/Disabled            |         |
| RS232 SELECTED     | Disabled      | Disabled                 | Enabled/Disabled            |         |
| NEOI WARNING       | 1min          | 5mins                    | 1min - 15mins               |         |
| EOI POINT          | 1.0%          | 1.0%                     | 0.1% - 5% of syringe volume |         |
| KVO AT EOI         | Enabled       | Enabled                  | Enabled/Disabled            |         |
| KVO RATE           | 1.0ml/h       | 1.0ml/h                  | 0.1ml/h - 2.5ml/h           |         |
| BACK OFF           | Disabled      | Enabled                  | Enabled/Disabled            |         |
| AC FAIL            | Enabled       | Enabled                  | Enabled/Disabled            |         |
| PRESSURE DISPLAY   | Disabled      | Enabled                  | Enabled/Disabled            |         |
| PRESSURE DEFAULT   | L-5           | L-3                      | L0 - 10(50mmHg -1000mmHg)   |         |
| WEIGHT             | 70.0Kg        | 70.0Kg                   | 0.01Kg - 250Kg              |         |
| PURGE RATE         | 200ml/h       | 200ml/h                  | 100ml/h - 500ml/h           |         |
| PURGE VOLUME LIMIT | 2.0ml         | 2.0ml                    | 0.5ml - 5.0ml               |         |
| PURGE SYRINGE      | Enabled       | Disabled                 | Enabled/Disabled            |         |
| HANDS FREE BOLUS   | Enabled       | Enabled                  | Enabled/Disabled            |         |
| DEFAULT BOLUS VOL  | 5.0ml         | 5.0ml                    | 0.1ml - 100ml               |         |
| DEFAULT BOLUS RATE | 1200ml/h      | 1200ml/h                 | 150ml/h - 1200ml/h          |         |
| MANUAL BOLUS       |               | Disabled                 | Enabled/Disabled            |         |
| CALLBACK TIME      |               | 2.0mins                  | 0.1mins - 15.0mins          |         |
| EVENT LOG DISPLAY  | Disabled      | Enabled                  | Enabled/Disabled            |         |
| BATTERY ICON       |               | Enabled                  | Enabled/Disabled            |         |
| AUDIO VOLUME       | Medium        | Medium                   | Low, Medium, High           |         |
| AUTO NIGHT MODE    | Enabled       | Enabled                  | Enabled/Disabled            |         |

### Syringes Enabled

| Make | Size(s) | Make | Size(s) |
|------|---------|------|---------|
|      |         |      |         |
|      |         |      |         |
|      |         |      |         |
|      |         |      |         |

**Hospital Name**

**Serial No.**

**Software Version**

**Approved by**

**Configured by**

**Date**

**Date**

## Alaris® CC Syringe Pump Drug Protocol Setup

---

---

[illegible]

---

---

---

---

---

\* - 100 drug names with a maximum of 17 characters are available for V2.3.x & above software.

## Alaris® Syringe Pumps

84/88

\* - 100 drug names with a maximum of 17 characters are available for V2.3.x & above software.

## 1000SM00001 Issue 9

# ***Appendix F***

## ***Service Contacts***

## Service Contacts

For service, contact your local Cardinal Health Affiliate Office or Distributor.

### AE

Cardinal Health,  
PO Box 5527,  
Dubai, United Arab Emirates.  
Tel: (971) 4 28 22 842  
Fax: (971) 4 28 22 914

### ES

Cardinal Health,  
Avenida Valdeparra 27,  
Edificio Alcor,  
28108 - Alcobendas, Madrid, España.  
Tel: (34) 91 657 20 31  
Fax: (34) 91 657 20 42

### NO

Cardinal Health,  
Solbråveien 10 A,  
1383 ASKER,  
Norge.  
Tel: (47) 66 98 76 00  
Fax: (47) 66 98 76 01

### AU

Cardinal Health,  
8/167 Prospect Highway,  
Seven Hills, NSW 2147,  
Australia.  
Tel: (61) 2 9838 0255  
Fax: (61) 2 9674 4444  
Fax: (61) 2 9624 9030

### FR

Cardinal Health,  
Immeuble Antares - Technoparc,  
2, rue Charles-Edouard Jeanneret.  
78300 POISSY,  
France.  
Tél: (33) 1 30 06 74 60  
Fax: (33) 1 39 11 48 34

### NZ

Cardinal Health,  
14 George Bourke Drive  
Mt Wellington, Auckland  
PO Box 14234  
Panmure, Auckland  
Tel: 09 270 2420  
Freephone: 0508 422734  
Fax: 09 270 6285

### BE

Cardinal Health,  
Otto De Mentockplein 19,  
1853 Strombeek - Bever,  
Belgium.  
Tel: (32) 2 267 38 99  
Fax: (32) 2 267 99 21

### GB

Cardinal Health,  
The Crescent, Jays Close,  
Basingstoke,  
Hampshire, RG22 4BS,  
United Kingdom.  
Tel: (44) 0800 917 8776  
Fax: (44) 1256 330 860

### SE

Cardinal Health,  
Hammarbacken 4B,  
191 46 Sollentuna,  
Sverige.  
Tel: (46) 8 544 43 200  
Fax: (46) 8 544 43 225

### CA

Cardinal Health,  
235 Shields Court,  
Markham,  
Ontario L3R 8V2,  
Canada.  
Tel: (1) 905-752-3333  
Tel: 800-908-9918  
Fax: (1) 905-752-3343

### HU

Cardinal Health,  
Döbrentei tér 1,  
H-1013 Budapest,  
Magyarország.  
Tel: (36) 14 88 0232  
Tel: (36) 14 88 0233  
Fax: (36) 12 01 5987

### US

Cardinal Health,  
10221 Wateridge Circle,  
San Diego, CA 92121,  
USA.  
Tel: (1) 800 854 7128  
Fax: (1) 858 458 6179

### CN

Cardinal Health,  
Shanghai Representative Office,  
Suite 9B, Century Ba-Shi Building,  
398 Huai Hai Rd(M.),  
Shanghai 200020,  
China.  
Tel: (56) 8621-63844603  
Tel: (56) 8621-63844493  
Fax: (56) 8621-6384-4025

### IT

Cardinal Health,  
Via Ticino 4,  
50019 Sesto Fiorentino,  
Firenze, Italia.  
Tel: (39) 055 30 33 93 00  
Fax: (39) 055 34 00 24

### ZA

Cardinal Health,  
Unit 2 Oude Molen Business Park,  
Oude Molen Road, Ndabeni,  
Cape Town 7405, South Africa.  
Tel: (27) 860 597 572  
Tel: (27) 21 510 7562  
Fax: (27) 21 5107567

### DE

Cardinal Health,  
Pascalstr. 2,  
52499 Baesweiler,  
Deutschland.  
Tel: (49) 2401 604 0  
Fax: (49) 2401 604 121

### NL

Cardinal Health,  
Kantoren pand "Hoefse Wing",  
Printerweg, 11,  
3821 AP Amersfoort,  
Nederland.  
Tel: (31) 33 455 51 00  
Fax: (31) 33 455 51 01

# ***Appendix G***

## ***Document History***

## Document History

| Issue | Date     | CO No. | Author    | Update Description                                                                                                                                                                                                                                                                                                                                                    |
|-------|----------|--------|-----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1     | 19/12/02 | 4091   | Ian Tyler | Initial release                                                                                                                                                                                                                                                                                                                                                       |
| 2     | 14/01/03 | 4268   | Ian Tyler | RS232 pin-out table pin descriptions pin5 and 9 switched.<br>Drug Protocol Setup table - Dose Rate Max & Min headings switched.                                                                                                                                                                                                                                       |
| 3     | 24/04/03 | 4432   | Ian Tyler | Note added on cleaning with reference to shelf keypads.<br>Added pole clamp pivot screw and changed one screw description.<br>Added pole clamp arm and updated part numbers.<br>Drawings of Pole clamp and Transducer to show latest versions.                                                                                                                        |
| 4     | 10/06/04 | 4710   | Ian Tyler | Information relating to Guardrails® Safety Software.<br>Pole clamp arm replacement.<br>Updates for latest Software versions, including - configuration options, part numbers, error codes and new screen display.<br>Enhanced 2 piece syringe flange clamp.<br>Adhesive feet availability information.<br>Added additional diagrams for Chassis and Plunger assembly. |
| 5     | 27/09/04 | 5478   | Ian Tyler | Administration change.                                                                                                                                                                                                                                                                                                                                                |
| 6     | 6/04/05  | 5688   | Ian Tyler | Add Alaris® PK Syringe Pump                                                                                                                                                                                                                                                                                                                                           |
| 7     | 20/04/05 | 5920   | Ian Tyler | Add more information regarding battery calibration<br>Updated TIVA Drug setup & protocol to include more units in column headings<br>Add access code 612                                                                                                                                                                                                              |
| 8     | 4/11/05  | 6064   | Ian Tyler | Update Speed test values<br>New PK error codes<br>Strain plate update<br>New software settings<br>Pole clamp part numbers update<br>Rebranded from ALARIS Medical Systems to Cardinal Health<br>Added access codes 175 & 711                                                                                                                                          |
| 9     | 23/06/06 | 6932   | Ian Tyler | Administration rebrand change.<br>Updated Motor Plate replacement instructions & drawings                                                                                                                                                                                                                                                                             |